

Avio ICP-OES product line



Helmut Ernstberger

13.6.2019, Budapest Seminar



PerkinElmer long history in ICP-OES

1979 – PerkinElmer Introduces the ICP 5000

1984 – PerkinElmer Plasma II Introduced

1987 – PerkinElmer introduces the Plasma40

1993 – PerkinElmer introduces Optima 3000

1995 – Optima 3000DV Introduced

2001 – Optima 4300 and Optima 2000

2005 – Optima 5300 and Optima 2100

2007 – Optima 7300 and Optima 7000

2011 – Optima 8x00 Family


2016 – Avio 200

2017 – Avio 500



PerkinElmer Continues to Lead the Market for ICP-OES

ICP-OES Vendor Share 2017



Source: SDI Global Assessment Report 2018

New AVIO series

- Avio 200, launched July 2016



Scanning array ICP-OES

- Avio 500, launched July 2017




Simultaneous ICP-OES

What do Avio 200 and 500 have in common?

- Flat Plate technology for low argon consumption
- 40 MHz RF generator for robustness and stability
- PlasmaShear
- Vertical torch
- Easy to use torch mount
- Dual View
- Viewing height selectable element by element
- PlasmaCam
- Interference correction capabilities IEC, MSF
- Syngistix Software
- RoHS compliant
- 4-channel, 12-roller peristaltic pump
- Compact size

Argon saving: Flat Plate Technology



Patented Flat Induction Plates operate at **40 – 60% of the Argon flow** of helical designs



- Avio and Optima platforms are the only ICP-OES systems capable of running a stable plasma at 8 L/min plasma gas
- Same robust plasma conditions for all samples
- Full power range allows the analysis of all sample types
- No helical load coil needed: eliminates maintenance, consumables cost, downtime and leakage risk

Load coil vs Flat plate – Performance evaluation

- Side by side comparison
 - Optima 8300 vs Optima 4300 DV
 - Difference: RF generator / Flat Plate
 - Same sample introduction system: GemCone / Cyclonic
- Matrix Robustness according to Mermet's criteria for Mg
 - Mg II (280.270 nm) / Mg I (285.213 nm) > 8 for robust plasma
 - Vary Neb flow and plasma power

Plasma parameters	Load Coil	Flat Plate
Plasma Power [W]	varied	varied
Plasma gas flow [L/min]	15	8
Aux. gas flow [L/min]	0.2	0.2
Neb. gas flow [L/min]	varied	varied


- High Matrix Samples
 - Digested coal fly ash: NIST 1633c
 - Mg in digest = 10 mg/L

Load coil vs Flat plate - Plasma robustness


UNIVERSITY OF JYVÄSKYLÄ

Radial measurement

a) load coil



b) flat plate




Determined Mg II / Mg I line intensity ratios (radial) in SRM1633c using a) load coil and b) flat plate plasma with GemCone nebulizer and cyclonic spray chamber

Better efficiency of energy transfer for Flat Plate


Load coil vs Flat plate - Sensitivity

UNIVERSITY OF JYVÄSKYLÄ

a) load coil



b) flat plate



Determined intensities (radial) of the Mg II 280.270 nm in SRM1633c using a) load coil and b) flat plate plasma

Higher sensitivity for Flat Plate

Load coil vs Flat plate – Method detection limits

- MDLs determined according to EPA 200.7
- Axial view
- Matrix: Aqua regia, 1:5 diluted

Element	Wavelength [nm]	Flat Plate MDL [µg/L]	Improvement factor over load coil
As	188.979	8.2	1.5
Co	228.616	1.6	19
Cr	267.716	1.5	23
Cu	324.752	2.7	10
Mn	257.610	2.7	7
Ni	231.604	0.6	23
Pb	220.353	2.1	9
Zn	213.857	0.5	44

Flat Plate detection limits on average 17 times lower

Source: University of Jyväskylä, with permission

Saving Argon: PlasmaShear



PlasmaShear exclusive to PerkinElmer

- Improved performance
- Efficient removal of the tail plume **without the need of argon**
- Eliminates self-absorption
- Extended linear range
- Matrix independent
- Increase accuracy
- Reduces maintenance time

Vertical Torch Configuration

New Vertical Torch Configuration

- Easy-to-use torch mount with fewer parts
 - No tools required
 - One-handed removal
 - Quick-change components
- Easy and cost-effective maintenance
 - Removable injector independent of torch
- Lockable depth position for repeatable torch insertion
- Adjustable, even while the ICP is running, to ensure optimum results
- Solid design
- Accommodates a vast array of sample introduction accessories that make running sample analysis simple, regardless of sample matrix




Dual View: the Best of Both Worlds – No Compromise


Designers of the first Dual View ICP-OES (Optima 3000 in 1995)

- Provides analysis of both low-level and high-level elements in one method
 - Sub-ppb to % level concentration range for all elements across UV and Vis
- Fully-selectable radial viewing height and axial sampling depth
- Delivers flexibility for different sensitivities; i.e. Na and K
- Allows regulatory methods to be run at prescribed wavelength
- Removes EIE interferences
- 1500W in both axial and radial modes

Lower Cost and Time Savings




PlasmaCam Viewing Camera



- Real-time plasma viewing
- Verify sample introduction components during analysis
- Simplifies method development
- Perform remote diagnostics
- Ensures maximum up time

Industry's first integrated color camera for continuous plasma viewing during analysis


Adjustable viewing height can eliminate EIE Effects



	F'n	Element	Plasma (L/min)	Aux (L/min)	Neb (L/min)	Power (watts)	View Dist	Plas View
10	A	K 766.490 rad I+Srl	15	0.2	0.70	1400	15.0	Radial
11	A	K 766.490 rad I+Y	15	0.2	0.70	1400	10.0	Radial
12	A	Mg 279.077 rad II+Y	15	0.2	0.70	1400	15.0	Radial
13	A	Mg 285.213 rad I+In	15	0.2	0.70	1400	15.0	Radial
14	A	Mn 260.568 ax II+Y	15	0.2	0.70	1400	15.0	Axial
15	A	Mn 257.610 ax II+Y	15	0.2	0.70	1400	15.0	Axial
16	A	Na 588.995 rad I+Srl	15	0.2	0.70	1400	10.0	Radial
17	A	Na 589.592 rad I+Srl	15	0.2	0.70	1400	15.0	Radial
18	A	P 213.617 rad I+In	15	0.2	0.70	1400	15.0	Radial
19	A	P 213.617 ax I+In	15	0.2	0.70	1400	15.0	Axial


Avio 500 ICP-OES Features

- Simultaneous instrument
- High performance optical system
- Universal Data Acquisition
- Argon Quick Change
- SmartRinse




Analysis Speed: Optics & Detector

- Unique optical system & detector allow simultaneous reading of both the visible and UV regions of the spectrum
- Requires thermally stable optical bench



Avio 500 Performance: Unsurpassed stability




Meinhard Type K1 and standard baffled cyclonic spray chamber

9-hour stability run results in <1% RSD

Universal Data Acquisition - UDA

- UDA gathers data for all 6300+ wavelengths available on **both** detectors for each sample
- Simple click to turn on
- Uses approx 80kB
- No time penalty
- Live concentration results on screen for every element



Benefits

- No need to re-analyse samples for additional elements/ to confirm results
- Saves analyses time (and money) during method development, all of this can be done offline (not using gas, sample.)

UDA for Quantitative Analysis Scenario...

Method Editor : Effluent analyses

Spectrometer | Sampler | Process | Calibration | Checks | QC | Options

Define elements

Method description

	Symbol	Wavelength (nm)	Name	Function
1	Al	396.153	Al 396.153	Analyte
2	As	188.979	As 188.979	Analyte
3	Be	313.107	Be 313.107	Analyte
4	Cd	228.802	Cd 228.802	Analyte
5	Cd	214.440	Cd 214.440	Analyte
6	Co	228.616	Co 228.616	Analyte
7	Co	231.160	Co 231.160	Analyte
8	Cr	267.716	Cr 267.716	Analyte
9	Cu	327.393	Cu 327.393	Analyte
10	Cu	324.752	Cu 324.752	Analyte
11	Fe	238.204	Fe 238.204	Analyte
12	Mn	257.610	Mn 257.610	Analyte
13	Mo	202.031	Mo 202.031	Analyte
14	Ni	231.604	Ni 231.604	Analyte
15	Pb	220.353	Pb 220.353	Analyte
16	Y	371.029	Y 371.029	Int. Std.
17	Zn	206.200	Zn 206.200	Analyte
18				
19				

Define Elements


Settings

Spectral Windows

Periodic table

Wavelength

Elements and wavelengths can be selected by clicking on one of the buttons to the right



Data Viewer

Corrected Intensities | Conc. in Calib. Units | Conc. in Sample Units | Internal Standards | QC

☐ Show RSDs

	Sample Id	Y 371.0...	Al 396.1... (mg/L)	As 188.... (mg/L)	Be 313.... (mg/L)	Cd 228.... (mg/L)	Cd 214.... (mg/L)	Co 228.... (mg/L)	Co 231.... (mg/L)	Cr 267.... (mg/L)	Cu 327.... (mg/L)	Cu (n
1	blank	100.0%	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000
2	Trade Std#1	99.1%	10.000	4.000	0.400	0.400	4.000	4.000	4.000	4.000	4.000	4.000
3	Trade Std#2	98.9%	30.000	12.000	1.200	1.200	12.000	12.000	12.000	12.000	12.000	12.000
4	Trade Std#3	99.2%	50.000	20.000	2.000	2.000	20.000	20.000	20.000	20.000	20.000	20.000
5	Trade AQC	98.0%	25.319	10.084	1.046	1.067	10.252	10.107	10.141	10.350	10.099	10.099
6	TW1 US	106.9%	16.313	-0.320	-0.002	-0.009	0.001	0.060	0.043	0.017	0.320	0.117
7	TW1 SP	106.0%	28.386	4.328	0.442	0.407	4.224	4.100	4.134	4.183	4.345	4.224

Method Editor : Effluent analyses

Spectrometer | Sampler | Process | Calibration | Checks | QC | Options

Define elements

Method description

	Symbol	Wavelength (nm)	Name	Function
1	Al	396.153	Al 396.153	Analyte
2	As	188.979	As 188.979	Analyte
3	As	193.696	As 193.696	Analyte
4	Be	313.107	Be 313.107	Analyte
5	Cd	214.440	Cd 214.440	Analyte
6	Cd	228.802	Cd 228.802	Analyte
7	Co	231.160	Co 231.160	Analyte
8	Co	228.616	Co 228.616	Analyte
9	Cr	267.716	Cr 267.716	Analyte
10	Cu	324.752	Cu 324.752	Analyte
11	Cu	327.393	Cu 327.393	Analyte
12	Fe	238.204	Fe 238.204	Analyte
13	Mn	257.610	Mn 257.610	Analyte
14	Mo	202.031	Mo 202.031	Analyte
15	Ni	231.604	Ni 231.604	Analyte
16	Pb	220.353	Pb 220.353	Analyte
17	Y	371.029	Y 371.029	Int. Std.
18	Zn	206.200	Zn 206.200	Analyte
19				

Elements and wavelengths can be selected by clicking on one of the buttons to the right

Periodic table

Wavelength table

Define Elements

Settings

Spectral Windows

Cr 267.... (mg/L)	Cu 327.... (mg/L)	Cu 324.... (mg/L)
0.000	0.000	0.000
4.000	4.000	4.000

Periodic Table

Enter in Method

λ Table

Active row in Method Editor 19

Element Arsenic (As)

Wavelength

188.979
193.696
197.197
228.812

Enter calibration concentrations and reprocess data

Method Editor : Effluent analyses

Spectrometer | Sampler | Process | Calibration | Checks | QC | Options

Calibration units and standard concentrations

	Analyte	Calib Units	Trade Std#1	Trade Std#2	Trade Std#3
1	Al 396.153	mg/L	10	30	50
2	As 188.979	mg/L	4	12	20
3	As 193.696	mg/L	4	12	20
4	Be 313.107	mg/L	0.4	1.2	2
5	Cd 214.440	mg/L	4	12	20
6	Cd 228.802	mg/L	0.4	1.2	2
7	Co 231.160	mg/L	4	12	20
8	Co 228.616	mg/L	4	12	20
9	Cr 267.716	mg/L	4	12	20
10	Cu 324.752	mg/L	4	12	20
11	Cu 327.393	mg/L	4	12	20
12	Fe 238.204	mg/L	10	30	50
13	Mn 257.610	mg/L	4	12	20
14	Mo 202.031	mg/L	4	12	20
15	Ni 231.604	mg/L	4	12	20
16	Pb 220.353	mg/L	4	12	20
18	Zn 206.200	mg/L	10	30	50

Define Standards
Calib Units and Concentrations
Blank Usage
Equations and Sample Units
Initial Calibration
Multiline Calibration


Quantitation with UDA standard sets

- UDA Standard set
 - Covers a wide element range (69 elements)
 - Include UDA standards in run
- ALL wavelengths info is collected with UDA
- Retrospective quant analysis
 - Calibrate any UDA element not part in your original standards mix
 - Assign internal standard
 - Reprocess with UDA solution as standard
 - Quantitative analysis result – however only 1 concentration level

Answer questions about your samples which are asked in hindsight!

Argon Quick-Change

- A mode in the software that allows Ar tanks to be changed without powering off the Avio 500
- Instrument goes into Stand-by Mode
- The user has 30 minutes to change the Ar tank
- After the tank is changed, start up time is 15 minutes
 - Time required to purge the detector



Analysis Speed

- Sample: 1 ppm mixed element standard
- Number of wavelengths: 20
 - **Elements:** Al, As, Ba, Be, Ca, Cd, Co, Cr, Cu, Fe, K, Mg, Mn, Na, Ni, Pb, Se, Sr, V, Zn
 - Some elements with low sensitivity (i.e. As, Se, Pb, etc.)
 - Some elements with high sensitivity (i.e. Mg, Na, Sr, etc.)
- Default Read Time Range: 1 - 5 seconds
 - Software automatically selects best integration time within the range based on a pre-shot
 - Time selected per wavelength
 - Integration Time based on signal intensity
- Views: axial & radial
- **Analysis Time: 30 seconds**

Analysis Speed: Sample Washout

- Sample washout can be the longest time going from sample-to-sample
- A variety of sample introduction components are available to minimize washout
 - ESI FAST
 - Cetac ASXpress
 - Glass Expansion Niagara
- A variety of rinse options are available in Syngistix

Sample Washout: SmartRinse

- User-defined concentrations for each element in the method
- When the concentrations are below the user-defined thresholds, rinsing is complete

The screenshot shows the 'Method Editor : k-2007-Normal' window with the 'Sampler' tab selected. The 'Autosampler wash' section has 'SmartRinse' selected. A 'SmartRinse Concentration Limits' dialog box is open, displaying a table of target concentrations for various analytes.

Method Editor : k-2007-Normal

Spectrometer | **Sampler** | Process | Calibration | Checks | QC | Options

Autosampler wash

☒ SmartRinse Set...

☐ Never

☐ Between samples

☐ Only after samples whose concentrations exceed limits Set...

☐ After every sample + extra time if sample concentrations exceed limits Set...

Rate: 1.00 mL/min

Wash location: 0

Max time - samples: 120 sec

Wash time - calib. solutions: 30 sec

SmartRinse Concentration Limits

Target concentrations to reach before ending wash

	F'n	Analyte	Units	Concentration
1	A	Ag 328.068	mg/L	
2	A	Al 396-rad	mg/L	
3	A	Al 308-rad	mg/L	
4	A	Al 394-rad	mg/L	
5	A	As 188.979	mg/L	
6	A	As 193.696	mg/L	
7	A	B 249.677	mg/L	
8	A	Ba 493.408	mg/L	
9	A	Be 313.107	mg/L	

Concentration units: ☒ Calibration ☐ Sample

OK Cancel

Sample Washout: Upper Limits-1

- User-defined upper-limit concentrations for each element to determine if washing is required

Method Editor : k-2007-Normal

Spectrometer | **Sampler** | Process | Calibration | Checks | QC | Options

Autosampler wash

☐ SmartRinse Set...

☐ Never

☐ Between samples

☒ Only after samples whose concentrations exceed limits Set...

☐ After every sample + extra time if sample concentrations exceed limits Set...

Rate: 1.00 mL/min

Wash location: 0

Normal time: 60 sec

Extra time: 60 sec


Extra Wash Concentration Limits

	Analyte	Calibration Units	Concentration
1	Ag 328.068	mg/L	5
2	Al 396-rad	mg/L	5
3	Al 308-rad	mg/L	0
4	Al 394-rad	mg/L	0
5	As 188.979	mg/L	5
6	As 193.696	mg/L	5
7	B 249.677	mg/L	5
8	Ba 493.408	mg/L	5
9	Be 313.107	mg/L	5
10	Ca 315-rad	mg/L	5

OK Cancel

Sample Washout: Upper Limits-2

- User-defined upper-limit concentrations for each element to determine if extra washing is required



Full Range of Intelligent Autosamplers and Diluters



prep3

- Automatic calibration from stock standard solution(s)
- Automatic sample dilution
- Automatic over-range dilution
- Automatic internal standard out-of-range dilution
- Analyze samples from undiluted to up to 25x diluted
- Fully supported by PerkinElmer Syngistix software



<https://www.youtube.com/watch?v=yvWly-fDyuA>

P1 = Sample, P2 = Diluent, P3 = Internal Standard

$P1 + P2 = \text{Constant}$, $P1/(P1+P2) = \text{Dilution Factor}$



AMAZINGLY
CAPABLE
REMARKABLY
AFFORDABLE



AVIO 200 ICP-OES

Avio 200...

- Scanning Array Spectrometer
- Does not require to be thermostated
- Designed to turn off when not in use
- OFF is OFF
- Optical system with very high light throughput allows short read times
- Fastest start up from cold (10 min)
- DWS for best stability
- Smallest ICP on the market at 65x76x81cm



Fastest Start up



- 10 minutes from cold
 - No gas
 - No electricity
- How? DWS – Dynamic Wavelength Stabilisation
 - Use of the Neon reference source corrects for any peak shift due to temperature and purge (or lack of)

Benefits


- Excellent stability despite changing laboratory environment
- Best precision and accuracy

Dynamic Wavelength Stabilisation (DWS)


- All wavelengths are automatically referenced to Neon to eliminate wavelength drift
- done in real time for every replicate of every measurement
- Temperature fluctuations allowed



With and Without DWS



With: 20 hour run over
a temperature range of
10-35 °C



Without: 20 hour run
run over a temperature
range of 10-35 °C

Avio 200 Speed – 6 wavelengths, radial

Sample to sample times with autosampler without FAST valve

	Read time 0.5-2s	Read time 0.2-0.2s
Flush and read delay (typical)	40-50 s	40-50 s
Analytical time	72 s	48
Washout (typical)	20 s	20 s
Sample to sample time	137	113

Sample to sample times with autosampler with FAST valve

	Read time 0.5-2s	Read time 0.2-0.2s
Flush and read delay (typical)	10-15 s	10-15 s
Analytical time	72 s	48
Washout (typical)	10-15 s	10-15 s
Sample to sample time	97	73


Analytical time is average for solutions with range 0 – 100 ppm analyte

Attenuation mode

- Ability to extend the dynamic range is valuable for the analysis of high concentration samples.
- Attenuation Mode provides a simple way to dramatically extend the dynamic range (ca. 50x)
- Eliminate the need for additional sample prep or repeat analyses.
- Avoids the efforts involved with validating alternate wavelengths.
- Freely selectable, element by element.
- Allows measurement of both high & low concentration elements in the same method, using a mix of attenuated and unattenuated modes.

Attenuation mode

- Comparison of detector saturation limits for Attenuated and Normal modes Axially and Radially



Attenuation Mode allows concentrations greater than 1000 mg/L to be measured without dilution

The Avio Series

- **When to choose an Avio 500?**
 - Simultaneous instrument
 - High sample throughput regardless of number of wavelengths
 - Instrument used continuously or in planned intervals - Continuously ON
 - Backup quantitation with UDA
- **When to choose an Avio 200?**
 - Scanning array instrument
 - High sample throughput for < 10 wavelengths
 - medium sample throughput for > 10 wavelengths
 - Instrument used intermittently or ad hoc – OFF is OFF, short startup
 - versatility



Thank you for your attention!

Helmut Ernstberger

Product Specialist Inorganic EMEA EMT and Spain
helmut.ernstberger@perkinelmer.com