



Instructions for PerkinElmer® New Coronavirus Nucleic Acid Detection Kit

v 1.0

For prescription use only. For in vitro diagnostic use only.
For Emergency Use Authorization only.

This manual is proprietary to PerkinElmer, Inc., and intended only for customer use in connection with the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose without the prior written consent of PerkinElmer. **Follow the protocol included with the kit.**

Table of Contents

Product Name	3
Intended Use	3
Principles of the Assay	3
Kit Components and Packaging Specifications	4
Materials Required but Not Provided	4
Storage & Handling Requirements	5
Warnings and Precautions	5
Instruments	6
Collection, Storage & Shipment of Specimens	6
Assay Procedure	7
Interpretation of Results	10
Kit Limitations	12
Conditions of Authorization for the Laboratory	13
Assay Performance	14
Limit of Detection.....	14
LoD Verification.....	15
Analytical Reactivity (Inclusivity)	16
Analytical Specificity (Cross-reactivity)	16
Clinical Study.....	19
References	21

Product Name

PerkinElmer® New Coronavirus Nucleic Acid Detection Kit

Kit Contents

48 Tests

Intended Use

The PerkinElmer® New Coronavirus Nucleic Acid Detection Kit test is a real-time RT-PCR *in vitro* diagnostic test intended for the qualitative detection of nucleic acid from the SARS-CoV-2 virus in human oropharyngeal swab and nasopharyngeal swab specimens collected from individuals suspected of COVID-19 by their healthcare provider. Testing is limited to laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, to perform high complexity tests.

Results are for the identification of SARS-CoV-2 RNA. SARS-CoV-2 RNA is generally detectable in human oropharyngeal swab and nasopharyngeal swab specimens during the acute phase of infection. Positive results are indicative of presence of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. Laboratories within the United States and its territories are required to report all positive results to the appropriate public health authorities.

Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.

The PerkinElmer® New Coronavirus Nucleic Acid Kit is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures. The PerkinElmer® New Coronavirus Nucleic Acid Kit is only for use under the Food and Drug Administration's Emergency Use Authorization.

Principles of the Assay

The PerkinElmer® New Coronavirus Nucleic Acid Detection kit uses TaqMan-based real-time PCR technique to conduct *in vitro* reverse transcription of SARS-CoV-2 RNA, DNA amplification and fluorescence detection.

The assay targets specific genomic regions of SARS-CoV-2: nucleocapsid (N) gene and ORF1ab.

The TaqMan probes for the two amplicons are labeled with FAM and ROX fluorescent dyes respectively to generate target-specific signal.

The assay includes an RNA internal control (IC, bacteriophage MS2) to monitor the processes from nucleic acid extraction to fluorescence detection. The IC probe is labeled with VIC fluorescent dye to differentiate its fluorescent signal from SARS-CoV-2 targets.

The assay also uses a dUTP/UNG carryover prevention system to avoid contamination of PCR products and subsequent false positive results.

Kit Components and Packaging Specifications

Catalog Number: 2019-nCoV-PCR-AUS (48 tests/kit)

Component Name	Specifications & Loading		Main Ingredients	Storage Conditions
nCoV reagent A	950 µL	× 1 tube	Buffers, dNTPs, Mg ²⁺	-25 to -15°C
nCoV reagent B	230 µL	× 1 tube	TE Buffer, primers, probes	-25 to -15°C
nCoV enzyme mix	170 µL	× 1 tube	Taq DNA polymerase, MMLV, RNasin, UNG	-25 to -15°C
nCoV internal control	390 µL	× 1 tube	TE Buffer, bacteriophage MS2	-25 to -15°C
nCoV positive control	1.4 mL	× 2 tubes	SARS-CoV-2 RNA fragments capsulated in bacteriophage	-25 to -15°C
nCoV negative control	1.4 mL	× 2 tubes	TE Buffer	-25 to -15°C

Notes: 1) The reference materials and other components in the kit should be treated as potential sources of infection. 2) The use of this kit should be strictly in accordance with the nucleic acid amplification guidelines to operate in compliance with the requirements of the appropriate laboratories. 3) The components in different batches of the kit cannot be used interchangeably.

Materials Required but Not Provided

1. RNA extraction reagents: The PerkinElmer® Nucleic Acid Extraction Kits (KN0212) and PreNAT II (SY61)(software version 1.00.06).
2. Instrument and software: Applied Biosystems™ 7500 Real-Time PCR System (4351104 with Laptop, 4351105 with desktop) (software version 2.3).
3. Additional tools and consumables required for automatic nucleic acid extraction and PCR setup using Pre-NAT II:

Items	Cat. No.
Centrifuge	TDL-80-2B
Vortex mixer	XW-80A
900 µL conductive tip Sterilized	AF01MP-9-XS
175 µL conductive tip Sterilized	AF200P-9-XS
50 µL conductive tip Sterilized	ATO5OP-9-XS-LB
150 mL Reagent Trough	C3040016
33 mL Reagent Trough	CJ222161115
2 mL U type 96 deep-well plate	DP20UR-9-N
Rod sleeve	CMG-550
1.5 mL transparent centrifugal tube	MCT-150-C

0.2ml PCR 8-trip tubes	PCR-0208-C
Caps for 0.2ml PCR 8-trip tubes	PCR-2CP-RT-C
Deep-well plate sealing film	HY3020011

Storage & Handling Requirements

1. Store all reagents at -25 to -15°C.
2. Use the reagents within 30 days once opened.
3. Completely thaw the reagents before use.
4. Avoid excessive freeze/thaw cycles for reagents.

Warnings and Precautions

1. For in vitro diagnostic use under Emergency Use Authorization only.
2. Positive results are indicative of the presence of SARS-CoV-2 RNA.
3. Laboratories within the United States and its territories are required to report all positive results to the appropriate public health authorities.
4. Keep the kit upright during storage and transportation.
5. Before using the kit, check tubes for leakage or damage. Each component in the kit should be thawed at room temperature, thoroughly mixed, and centrifuged before use.
6. Cross-contamination may occur when inappropriate handling of reference materials and specimens, which will cause inaccurate results. It is recommended to use sterile disposable filter-tips to aspirate reagents and specimens.
7. All specimen to be tested and the reference materials of the kits should be considered as infectious substances and processed strictly in accordance with laboratory biosafety requirements. Sterile centrifuge tubes and filter-tips should be used. After use, the tips should be disposed into a waste bin containing a 10% sodium hypochlorite solution. After the operation, the work area surface and the instrument surface should be disinfected with a freshly prepared 10% sodium hypochlorite solution, and then cleaned with 75% ethanol or pure water. Finally, turn on UV light to disinfect working surfaces for 30 minutes.
8. The PCR instrument used for this assay should be calibrated regularly according to instrument's instructions to eliminate cross-talks between channels.
9. This kit uses PCR-based technology and experiments should be conducted in three separate areas: reagent preparation area, specimen preparation area, amplification area. Access to each area must be in strict accordance with a single flow direction, namely the specimen preparation area → reagent preparation area → amplification area. Protective equipment accessories (goggles, work clothes, hats, shoes, gloves, etc.) should be worn during operation

and protective equipment accessories should be changed when entering and leaving different work areas. Protective equipment accessories in each work area are not interchangeable.

Instruments

PerkinElmer® PreNAT II Automated Workstation
Applied Biosystems® 7500 Real-Time PCR system.

Collection, Storage & Shipment of Specimens

1. Respiratory Specimens

Use only synthetic fiber swabs with plastic shafts. Do not use calcium alginate swabs or swabs with wooden shafts, as they may contain substances that inactivate some viruses and inhibit PCR testing. Place swabs immediately into sterile tubes containing 3 ml of viral transport media. For initial testing, nasopharyngeal swab specimens are recommended. Collection of oropharyngeal swabs is a lower priority and is acceptable if other swabs are not available.

1.1 Nasopharyngeal swab (NP): Insert a swab into nostril parallel to the palate. Swab should reach depth equal to distance from nostrils to outer opening of the ear. Leave swab in place for several seconds to absorb secretions. Slowly remove swab while rotating it.

1.2 Oropharyngeal swab (e.g., throat swab, OP): Swab the posterior pharynx, avoiding the tongue.

2. Storage

Store specimens at 2-8°C for up to 72 hours after collection. If a delay in testing or shipping is expected, store specimens at -70°C or below.

3. Shipping

Specimens PUI's must be packaged, shipped, and transported according to the current edition of the International Air Transport Association (IATA) Dangerous Goods Regulation External Icon. Store specimens at 2-8°C and ship overnight to the lab on ice pack. If a specimen is frozen at -70°C ship overnight to the lab on dry ice. Additional useful and detailed information on packing, shipping, and transporting specimens can be found at Interim Laboratory Biosafety Guidelines for Handling and Processing Specimens Associated with Coronavirus Disease 2019 (COVID-19).

4. For more information, refer to:

Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens from Persons for Coronavirus Disease 2019 (COVID-19)

<https://www.cdc.gov/coronavirus/2019-nCoV/guidelines-clinical-specimens.html>

5. Interim Laboratory Biosafety Guidelines for Handling and Processing Specimens Associated with Coronavirus Disease 2019 (COVID-19)

<https://www.cdc.gov/coronavirus/2019-nCoV/lab-biosafety-guidelines.html>

Assay Procedure

Nucleic acid extraction and PCR set up using Pre-NAT II Automated Workstation

Pre-NAT II Automated Workstation is designed to process 1-96 samples for downstream molecular assays. It contains a liquid handling system which automatically pipettes and mixes reagents and samples, a purification module that extracts and purifies nucleic acids, and an automatic PCR setup function which is also conducted by the liquid handling system. The entire workflow is automatic without manual intervention. Detailed operation instructions of Pre-NAT II can be found in the Pre-NAT II Automated Workstation User Manual. A quick-start instruction for the SARS-CoV-2 Real-time RT-PCR assay is described as below.

1. Take the nCoV Internal Control, nCoV Positive Control and nCoV Negative Control out from freezer, place them in a biological safety cabinet and completely thaw them at room temperature. Vortex the tubes to mix the contents, then centrifuge the tubes briefly at 1000 rpm to collect the liquid to the bottom of the tubes.
2. Prepare specimens and place them in a biological safety cabinet. If the specimens are frozen, completely thaw them at room temperatures and follow the operations described in 1.1 for the controls.
3. Take the Magnetic Beads from the PerkinElmer® Nucleic Acid Extraction Kits (KN0212) kit, vortex the tube for one minute to completely suspend the beads in the solution.
4. Turn on the PreNAT II instrument, double click the “Pre-NAT II” software icon, select username and enter password to start, then follow software guidance to initialize the instrument.
5. After initialization, click “Program Input” to choose an extraction protocol. For the SARS-CoV-2 assay, choose “2019-nCoV” from the protocol list.
6. In the same window, input the number of specimens that are going to be processed at the indicated box, positive control and negative control should not be counted, as they are pre-set in the 2019-nCoV protocol. After the sample number is entered, click “Set Complete” to proceed to the loading guidance for reagents and consumables.
7. Remove the lids from reagents, controls and specimens, load the

consumables, reagents, specimens, and controls according to software guidance, then double check to confirm that all items are at the positions indicated by software. Close instrument door after finish loading. Click “run” to start the protocol, the procedures automatically performed by Pre-NAT II are described below.

- 7.1 Add 400 µL of each specimen, Negative Control and Positive Control to the wells of a 96 deep-well plate, and add 5 µL nCoV Internal Control, Lysis/Binding Buffer and Magnetic Beads to each well.
- 7.2 Magnetic rods take rod sleeves and rotate in 96 deep-well plate to mix (magnetic force off status), during which stage DNA/RNA is released through lysis and binds to magnetic beads.
- 7.3 During lysis and binding, automatic liquid handler pipettes Wash Buffer A to a 96 deep-well plate.
- 7.4 Magnetic force is turned on for magnetic rods and beads are collected from Lysis/Binding reaction to Wash Buffer A.
- 7.5 Magnetic rods (magnetic force off) rotate to wash beads in Wash Buffer A and proceed in a same manner to wash beads in Wash Buffer B.
- 7.6 Finally, the beads are collected and placed into elution buffer to elute DNA/RNA.
- 7.7 During elution, liquid handler pipettes/mixes PCR reagents to prepare a PCR mix and aliquot to PCR tubes. Elution volume is 40 µL for each sample.
- 7.8 Eluted DNA/RNA is added to PCR mix in each tube, which is ready for amplification.

If PCR setup needs to be done manually, please refer to the following procedures.

1. Prepare PCR mix in Reagent Preparation Area according to the following table. It is recommended to prepare 110% of the calculated amount of PCR mix to account for pipetting carryovers.

Component	Volume/ test	Volume for N Samples and 2 Controls	110% of volume
nCoV Reagent A	15 µL	15 x (n + 2) µL	16.5 x (n + 2) µL
nCoV Reagent B	3 µL	3 x (n + 2) µL	3.3 x (n + 2) µL
nCoV Enzyme mix	2 µL	2 x (n + 2) µL	2.2 x (n + 2) µL

2. Completely vortex the prepared PCR mix, aliquot 20 μL into each PCR tube or each well of a 96-well PCR plate.
3. Add 40 μL of extracted nucleic acid into each tube or well containing PCR mix, close lids for the PCR tubes or seal PCR plates with an appropriate film, slightly vortex the tubes and briefly centrifuge them to get rid of bubbles.

Amplification (in PCR area)

1. Set up and run the Applied Biosystems™ 7500 Real-Time PCR instrument. Refer to Applied Biosystems™ 7500 Real-Time PCR Instrument Reference Guide for detailed instructions. In general, double-click 7500 software 2.3 › New experiments › Setup Experiment Properties › Setup the Targets and Samples in Plate Setup › Setup Run Method, then click Run and Start.
2. When setup Experiment Properties, please check the following run settings and choose the correct settings.
 - Instrument: 7500 (96 wells)
 - Run type: Quantitation - Standard Curve
 - Run reagent: TaqMan reagents
 - Run mode: Standard
3. When setting up the Targets and Samples, create the following detectors with the quencher set as none. The passive reference must be set as None.

Target Name or Detector	Reporter	Quencher
N	FAM	None
ORF1ab	ROX	None
IC	VIC/HEX	None

4. Set up the plate layout by assigning a unique sample name to each well.
5. Assign a Task to each well.
 - Unknown: for patient samples
 - Standard: for Positive Control
 - NTC: for Negative Control
6. Set Run method as following for PCR amplification and fluorescence detection, the sample volume is 60 μL .

Step	Temperature	Time	Number of Cycles
1	37°C	2 minutes	1
2	50°C	5 minutes	1
3	42°C	35 minutes	1
4	94°C	10 minutes	1
5	94°C	10 seconds	45
	55°C	15 seconds	
	65°C*	45 seconds	

* Collect fluorescence signal during the final 65°C step.

7. Double check all settings then click Run and Start to initialize amplification.

Interpretation of Results

1. Baseline and threshold setting


After the run completion, save and analyze the data according to PCR instrument instructions.

1.1 Set baseline for each target

View the baseline values, in the Graph Type drop-down list, select Linear. Select the Baseline check box to show the start cycle and end cycle. The horizontal part of the baseline is used for the baseline range, which normally starts from 3-5 cycles and ends at 15-20 cycles. Baseline setting is normally automatically done by instrument. It can also be manually adjusted to choose the horizontal part of the curve.

1.2 Set threshold for each target

View the threshold values, In the Graph Type drop-down list, select Linear. In the Target drop-down list, select N or ORF1ab or IC. Select the Threshold check box to show the threshold. Thresholds should be adjusted to fall within exponential phase of the fluorescence curves and above any background signal (refer to Figure A, B). The threshold value for different instruments varies due to different signal intensities.



- 1.3 Perform data analysis by clicking “Analyze” button of the software.
- 1.4 Output the data to csv file by the “export” function of the software.
- 1.5 Interpret the results based on the tables listed in “Quality Control” and “Examination and Interpretation of Patient Specimen Results”

2. Quality Control

The product provides negative control, positive control, and internal control to monitor the reliability of the results for the entire batch of specimens from sample extraction to PCR amplification. All test controls should be examined prior to interpretation of patient results. If the controls are not valid, the patient results cannot be interpreted.

- 2.1 Negative Control: both ORF1ab and N of SARS-CoV-2 must be not detected, and the Ct value of internal control should be ≤ 40 ;
- 2.2 Positive Control: both ORF1ab and N of SARS-CoV-2 must be detected and their Ct values should be ≤ 35 , the Ct value of internal control does not have to be ≤ 40 for positive control.
- 2.3 If negative and positive control results are not as described above, the test results of the entire batch are invalid.
- 2.4 Internal Control: If the result for a specimen is SARS-CoV-2 RNA not detected, the Ct value of the internal control must be ≤ 40 , otherwise the result of that specimen is inconclusive; if the result for a specimen is SARS-CoV-2 RNA detected, the Ct value of the internal control is not required to be considered valid. Positive control, negative control and IC in positive and negative control should meet the requirements listed in the below table to ensure valid results.

Control type	Ct		
	N (FAM)	ORF1ab (ROX)	IC (HEX/VIC)
Negative	Undetermined or > 42	Undetermined or > 42	Ct ≤ 40
Positive	≤ 35	≤ 35	/

/: No requirements on the Ct value.

3. Examination and Interpretation of Patient Specimen Results:

Assessment of clinical specimen test results should be performed after the positive and negative controls have been examined and confirmed to be valid and acceptable. If the controls are not valid, the patient results cannot be interpreted.

The table below lists the expected results for the kit with valid positive control and negative control:

Ct		Result interpretation
IC (VIC/HEX)	N(FAM), ORF1ab (ROX)	
≤40	Both targets Undet or >42	SARS-CoV-2 not detected
/	Both targets ≤ 42	SARS-CoV-2 detected
/	One of the targets ≤ 42	SARS-CoV-2 detected
>40 or Undet	Both targets Undet or >42	Invalid result, specimen needs to be re-tested from re-extraction or re-collected from patient for test.

/: No requirements on the Ct value. Undet: Undetermi

Kit Limitations

1. The use of this assay as an in vitro diagnostic under FDA Emergency Use Authorization (EUA) is limited to laboratories that are certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. § 263a, to perform high complexity tests.
2. This kit is used for qualitative detection of SARS-CoV-2 RNA from human oropharyngeal swab and nasopharyngeal swab. The results cannot directly reflect the viral load in the original specimens.
3. The PerkinElmer New Coronavirus Nucleic Acid Detection Kit performance has only been established with the specimen types described in the Intended Use section. Testing other types of specimen may cause inaccurate results.
4. The specimens to be tested shall be collected, processed, stored and transported in accordance with the conditions specified in the instructions. Inappropriate specimen preparation and operation may lead to inaccurate results.
5. Extraction and amplification of nucleic acid from clinical samples must be performed according to the specified methods listed in this procedure. Other extraction approaches and processing systems have not been evaluated.
6. Amplification and detection of SARS-CoV-2 with the PerkinElmer New Coronavirus Nucleic Acid Detection Kit has only been validated with the Applied Biosystems® 7500 Real-Time PCR instrument. Use of other instrument systems may cause inaccurate results.
7. The limit of detection (LoD) is determined based on a 95% confidence of detection. When SARS-CoV-2 presents at or above the LoD concentration in the test specimen, there will be a low probability that SARS-CoV-2 is not detected. When SARS-CoV-2 presents below the

LoD concentration in the test specimen, there will also be certain probability that SARS- CoV-2 can be detected.

8. Primers and probes for this kit target highly conserved regions within the genome of SARS- CoV-2. Mutations occurred in these highly conserved regions (although rare) may result in RNA being undetectable.
9. This kit uses an UNG/dUTP PCR products carryover prevention system which can prevent contamination caused by PCR products. However, in the actual operation process, the amplicon contamination can be avoided only by strictly following the instructions of PCR laboratories.
10. Negative results do not preclude SARS-CoV-2 infections and should not be used as the sole basis for treatment or other management decisions.
11. The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs have not been evaluated.
12. Laboratories are required to report all positive results to the appropriate public health authorities.

Conditions of Authorization for the Laboratory

The PerkinElmer New Coronavirus Nucleic Acid Detection Kit Letter of Authorization, along with the authorized Fact Sheet for Healthcare Providers, the authorized Fact Sheet for Patients, and authorized labeling are available on the FDA website: <https://www.fda.gov/MedicalDevices/Safety/EmergencySituations/ucm161496.htm>.

However, to assist clinical laboratories using the PerkinElmer New Coronavirus Nucleic Acid Detection Kit (“your product” in the conditions below), the relevant Conditions of Authorization are listed below:

- A. Authorized laboratories¹ using your product will include with result reports of your product, all authorized Fact Sheets. Under exigent circumstances, other appropriate methods for disseminating these Fact Sheets may be used, which may include mass media.
- B. Authorized laboratories using your product will use your product as outlined in the Instructions for Use. Deviations from the authorized procedures, including the authorized instruments, authorized extraction methods, authorized clinical specimen types, authorized control materials, authorized other ancillary reagents and authorized materials required to use your product are not permitted.
- C. Authorized laboratories that receive your product will notify the relevant public health authorities of their intent to run your product prior to initiating testing.
- D. Authorized laboratories using your product will have a process in place for reporting test results to healthcare providers and relevant public health authorities, as appropriate.
- E. Authorized laboratories will collect information on the performance of your product and report to DMD/OHT7-OIR/OPEQ/CDRH (via email: CDRH-

EUA-Reporting@fda.hhs.gov) and You (via email: COVID-19.TechnicalSupport@PerkinElmer.com) any suspected occurrence of false positive or false negative results and significant deviations from the established performance characteristics of your product of which they become aware.

- F. All laboratory personnel using your product must be appropriately trained in RT-PCR techniques and use appropriate laboratory and personal protective equipment when handling this kit, and use your product in accordance with the authorized labeling.

¹ The letter of authorization refers to, “United States (U. S.) laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, to perform high complexity tests” as “authorized laboratories.”

Assay Performance

Limit of Detection

Limit of detection (LoD) was determined as the lowest concentration of SARS-CoV-2 that at which the PerkinElmer® New Coronavirus Nucleic Acid Detection Kit can detect at a ≥95% positive rate. Samples were prepared using pooled clinical oropharyngeal swab specimen matrix collected from 12 individuals at 4 different time points giving. The pooled oropharyngeal swab matrix was tested using PerkinElmer® New Coronavirus Nucleic Acid Detection Kit and confirmed to be negative. In the first part of the study, a total of six 10-fold dilutions of known concentrations of inactivated SARS-CoV-2 virus (Isolate 2/231/human/2020/CHN) were prepared in negative clinical matrix and processed using the PerkinElmer nucleic acid extraction kit on the PreNAT II automated nucleic acid extraction system. Four PCR replicates per concentration were tested.

The results are summarized in the following tables.

Table: LoD study results from 10-fold dilution of virus stock.

Dilution Fold	N		ORF1a b		Mean Ct		
	Conc. (copies/ml)	Detection rate	Conc. (copies/ml)	Detection rate	N	ORF1ab	IC
1.0E+04	274	4/4	83.7	4/4	34.88	34.29	33.17
1.0E+05	27.4	4/4	8.37	3/4	38.74	37.67	33.27
1.0E+06	2.74	2/4	0.837	2/4	39.57	38.71	33.11
1.0E+07	0.274	1/4	0.0837	1/4	40.11	38.75	33.44
1.0E+08	0.0274	0/4	0.00837	0/4	/	/	32.68
1.0E+09	0.00274	0/4	0.000837	0/4	/	/	33.02
Negative	0	0/4	0.00	0/4	/	/	32.83

Based on the previous results, an additional eight 2-fold dilutions of known

concentrations of genomic RNA were prepared in negative clinical matrix. Twenty individual extraction replicates per dilution were tested. . The results are summarized in the following table.

Table: LoD study results from 2-fold dilution of virus stock.

Dilution Fold	N		ORF1a _b		Mean Ct		
	Conc. (copies/ml)	Detection rate	Conc. (copies/ml)	Detection rate	N	ORF1ab	IC
1.0E+04	274	20/20	83.7	20/20	34.95	35.48	31.55
2.0E+04	137	20/20	41.85	20/20	35.93	36.23	31.65
4.0E+04	68.5	20/20	20.93	20/20	36.91	37.10	31.70
8.0E+04	34.25	19/20	10.46	19/20	38.15	38.64	31.61
1.6E+05	17.13	18/20	5.23	13/20	38.80	39.48	31.60
3.2E+05	8.56	11/20	2.62	11/20	39.44	39.93	31.28
6.4E+05	4.28	8/20	1.31	7/20	40.26	40.44	31.41
1.28E+06	2.14	5/20	0.65	3/20	40.10	40.65	31.16
negative	0	0/20	0	0/20	/	/	31.15

The probit predicted 95% detection rate is presented in the below table.

Table: Probit predicted 95% detection rate using inactivated cultured SARS-CoV-2 (GenBank: MT135042.1).

Probit predicted 95% detection rate (copies/mL)	
N	ORF1ab
24.884 (95% CI: 17.032 – 57.917)	9.307 (95% CI: 7.428 – 13.003)

The probit-predicted LoD is estimated to result in approximately 3 copies/ PCR reaction.

LoD Verification

The LoD predicted by probit analysis was further verified by testing 20 extraction replicates of sample at 1xLoD concentration. Inactivated virus (GenBank: MT135042.1) was diluted into oropharyngeal swab matrix according to target ORF1ab LoD. 20 replicates of the sample were extracted using Nucleic Acid Extraction Kit (KN0212) on Pre-NAT II and tested using the PerkinElmer® New Coronavirus Nucleic Acid Detection Kit. The results are summarized in the following table. The probit predicted LoD was verified as the detection rate of N is 100% and ORF1ab is 95%.

Table: LoD verification results.

Conc.	Detection rate		Mean Ct		
	N	ORF1ab	N	ORF1ab	IC
1xLoD	100% (20/20)	95% (19/20)	38.39	38.11	31.18

Analytical Reactivity (Inclusivity)

BLASTn analysis queries alignments were performed with the SARS-CoV-2 ORF1ab and N oligonucleotide primer and probe sequences with all publicly available nucleic acid sequences for 2019-nCoV in GenBank to demonstrate the predicted inclusivity of the PerkinElmer New Coronavirus Nucleic Acid Detection Kit.

All the alignments show 100% identity to the available 2019-nCoV sequences.

Analytical Specificity (Cross-reactivity)

Cross-reactivity of the PerkinElmer New Coronavirus Nucleic Acid Detection Kit was evaluated using both *in silico* analysis and wet testing against normal and pathogenic organisms found in the respiratory tract.

BLASTn analysis queries of the PerkinElmer New Coronavirus Nucleic Acid Detection Kit primers and probes were performed against public domain nucleotide sequences with default settings. The database search parameters were as follows:

- The match and mismatch scores were 1 and -3, respectively.
- The penalty to create and extend a gap in an alignment was 5 and 2, respectively.
- The search parameters automatically adjusted for short input sequences and the expected threshold was 1000.

In summary no organisms, including other related SARS-coronaviruses, exhibited >80% homology to the forward primer, reverse primer, and probe for either the ORF1ab or N target. The results of the *in silico* analysis suggest the PerkinElmer New Coronavirus Nucleic Acid Detection kit is designed for the specific detection of SARS-CoV-2, with no expected cross reactivity to the human genome, other coronaviruses, or human microflora that would predict potential false positive RT-PCR results.

Wet testing against normal and pathogenic organisms of the respiratory tract was performed to confirm the results of the *in silico* analysis. Each organism identified in the table below was tested in triplicate with the PerkinElmer New Coronavirus Nucleic Acid Detection kit at the concentrations indicated. Each replicate was tested with a different reagent lot. All results were negative.

Table 5: Organisms tested for cross-reactivity with the PerkinElmer New Coronavirus Nucleic Acid Detection Kit

Pathogen	Source	Concentration	
		Evaluation	Unit
Human coronavirus 229E	ATCC VR-740™	2.8×10^2	TCID ₅₀ /mL
Human coronavirus OC43	ATCC VR-1558™	2.8×10^3	TCID ₅₀ /mL
Adenovirus type 3	ATCC VR-847™	$5.0 \times 10^{5.5}$	TCID ₅₀ /mL
Adenovirus type 2	ATCC VR-846™	5.6×10^4	TCID ₅₀ /mL
Adenovirus type 31	ATCC VR-1109™	1.6×10^6	TCID ₅₀ /mL
Adenovirus type 37	ATCC VR-929™	1.8×10^4	TCID ₅₀ /mL
Adenovirus type 51	ATCC VR-1603™	2.3×10^6	TCID ₅₀ /mL
Parainfluenza virus type 1	ATCC VR-94™	2.8×10^4	TCID ₅₀ /mL
Parainfluenza virus type 2	ATCC VR-92D™	0.303	ng/μL
Parainfluenza virus type 3	ATCC VR-93™	$5.0 \times 10^{4.5}$	TCID ₅₀ /mL
Parainfluenza virus type 4a	ATCC VR-1378™	2.8×10^4	TCID ₅₀ /mL
Parainfluenza virus type 4b	ATCC VR-1377™	1.6×10^3	TCID ₅₀ /mL
Influenza A virus (H1N1pdm09)	ATCC VR-1736™	2.6×10^3	PFU/mL
Influenza A virus (seasonal H1N1)	ATCC VR-1520™	$5.0 \times 10^{4.5}$	TCID ₅₀ /mL
Influenza A virus (H3N2)	ATCC VR-1679™	$5.0 \times 10^{3.5}$	TCID ₅₀ /mL
Influenza B virus	ATCC VR-1807™	7.6×10^2	PFU/mL
Enterovirus A71	ATCC VR-1432™	$5.0 \times 10^{5.5}$	TCID ₅₀ /mL
Enterovirus D68	ATCC VR-1823™	1.6×10^6	TCID ₅₀ /mL
Respiratory syncytial virus	ATCC VR-1400™	$5.0 \times 10^{3.5}$	TCID ₅₀ /mL
Rhinovirus B17	ATCC VR-1663™	2.0×10^6	PFU/mL
Rhinovirus A2	ATCC VR-482™	8.9×10^4	TCID ₅₀ /mL
<i>Chlamydia pneumoniae</i>	ATCC 53592™	2.9×10^5	IFU/mL
<i>Haemophilus influenzae</i>	ATCC 51907D™	10	μg/mL
<i>Streptococcus pyogenes</i>	ATCC 700294D-5™	7	μg/ml

<i>Streptococcus salivarius</i>	ATCC BAA-250D-5™	5.2	µg/ml
<i>Bordetella pertussis</i>	ZeptoMetrix Panel	Unknown	
Measles virus	National Standard for Influenza A/B Viral Nucleic Acids Detection Kit	Unknown	
Mumps virus		Unknown	
<i>Staphylococcus aureus</i>		Unknown	
Influenza A virus (H7N9)		Unknown	
<i>Mycoplasma pneumoniae</i>	ATCC 15531™	3.5 x 10 ⁶	
Human cytomegalovirus	Symbio	Unknown	
Hepatitis A virus	Symbio	1.84E+05	copies/mL
Hepatitis B virus	WHO NIBSC 10/266	9.55E+05	IU/mL
Hepatitis C virus	WHO NIBSC 14/150	1.00E+05	IU/mL
Human immunodeficiency virus type I (HIV-1)	WHO NIBSC 16/194	1.26E+05	IU/mL
Human immunodeficiency virus type II (HIV-2)	WHO NIBSC 08/150	1.00E+03	IU/mL
Epstein-barr virus	Symbio	1.46E+05	copies/mL
Cytomegalovirus	Symbio	1.15E+04	copies/mL

Interfering Substances Studies:

The potential interference of the substances listed below were tested in both the presence and absence of SARS-CoV-2 RNA with the PerkinElmer New Coronavirus Nucleic Acid Detection Kit. SARS-CoV-2 positive samples were prepared by mixing each of the potentially interfering substances with the assay positive control (synthetic SARS-CoV-2 ORF1ab and N RNA template encapsulated in MS2 bacteriophage) at approximately 3X the LoD. All positive and negative samples yielded expected results.

Table 6: Substances Tested for Interference with the PerkinElmer New Coronavirus Nucleic Acid Detection Kit

Substance	Concentration Tested	Substance	Concentration Tested
Valacyclovir	3.6 mg/mL	Saline	1 mg/mL
Entecavir	24.6 ng/mL	Beclomethasone dipropionate	22.5 µg/mL
Adefovir	90 ng/mL	Dexamethasone acetate	375 µg/mL
Ribavirin	5 mg/mL	Triamcinolone tablets	25 µg/mL
Acyclovir	3.6 mg/mL	Mometasone furoate	41.7 µg/mL
Azithromycin	1.35 mg/mL	Fluticasone propionate	1 mg/mL
Clarithromycin	30 µg/mL	Oxymetazoline	15% v/v

		hydrochloride	
Ciprofloxacin	7.5 µg/mL	Sulfur ointment	0.05% v/v
Telbivudine	15 µg/mL	Pharyngitis lozenges	0.05% v/v
Efavirenz	12.2 µg/mL	Chlorhexidine benzocaine	1.25 mg/mL
Tenofovir	1335 ng/mL	Menthol	5% v/v
Zanamivir	5 mg/mL	Rheumatoid factor	/
Mupirocin	0.02% w/v	Systemic Lupus Erythematosus	/
Tobramycin	0.6 mg/mL	Antinuclear antibody	/
Flunisolide	20 mg/mL	Hemoglobin	5 mg/mL
Budesonide	16.7 µg/mL	Human serum albumin	60 mg/mL
Bilirubin	0.6 mg/mL	Triglycerides	25 mg/mL
		Human genomic DNA	3 mg/mL

Clinical Study

The performance of The PerkinElmer® New Coronavirus Nucleic Acid Detection Kit was evaluated using contrived clinical oropharyngeal swabs and nasopharyngeal swabs. In total, 141 healthy individuals with no COVID-19 infection history, no COVID-19 symptoms and no contact with SARS-CoV-2 infected patients within in 14 days were recruited for the study. Both oropharyngeal swabs and nasopharyngeal swabs were collected from the 141 healthy individuals by trained personnel. Samples were immediately screened with The PerkinElmer® New Coronavirus Nucleic Acid Detection Kit and stored frozen until use.

The inactivated cultured virus (GenBank: MT135042.1) was spiked into 47 of the oropharyngeal swabs and 47 of the nasopharyngeal swabs at various (2xLoD, 4xLoD, 10xLoD, 20xLoD, 50xLoD, 100xLoD, 200xLoD, 250xLoD and 500xLoD, according to the LoD of target ORF1ab) Of the 47 contrived positive samples, 20 were spiked at concentrations equivalent to 2X the LoD, 20 were spiked with concentrations equivalent to 4X the LoD, and 7 were spiked with concentrations ranging from 10X LoD to 500X LoD. The remaining 94 oropharyngeal swabs and 94 nasopharyngeal swabs were tested as negative clinical samples.

The 141 oropharyngeal samples and 141 nasopharyngeal samples were tested in a blinded fashion (samples were prepared and capped, then all the tubes were mixed in a box and extracted using Nucleic Acid Extraction Kit (KN0212) and automatic extraction instrument Pre-NAT II in a random order). Testing was performed in a total of four RT-PCR runs with one positive and one negative control included per run. Results of the study are summarized below.

Table: Positive and negative control results from clinical evaluation.

Run	Control	N Ct	ORF1a b Ct	IC Ct	Pass
Run 1	Positive control	30.34	29.08	29.43	Yes
	Negative control	Undetermined	Undetermined	32.59	Yes
Run 2	Positive control	31.59	30.99	32.47	Yes
	Negative control	Undetermined	Undetermined	33.97	Yes
Run 3	Positive control	31.79	30.98	31.27	Yes
	Negative control	Undetermined	Undetermined	31.35	Yes
Run 4	Positive control	30.92	30.28	28.72	Yes
	Negative control	Undetermined	Undetermined	34.54	Yes

Table: Clinical evaluation with oropharyngeal samples.

SARS-CoV-2 concentration	Samples (N)	Detection rate		Mean Ct		
		N	ORF1ab	N	ORF1ab	IC
2xLoD	20	20/20	20/20	37.05	37.03	31.90
4xLoD	20	20/20	20/20	35.48	35.56	32.58
10xLoD	1	1/1	1/1	34.93	35.58	33.98
20xLoD	1	1/1	1/1	34.94	34.38	30.72
50xLoD	1	1/1	1/1	34.53	34.17	34.44
100xLoD	1	1/1	1/1	32.17	31.48	31.33
200xLoD	1	1/1	1/1	33.38	32.33	34.94
250xLoD	1	1/1	1/1	32.15	31.44	34.73
500xLoD	1	1/1	1/1	30.32	30.27	33.38
Negative	94*	0/94	0/94	/	/	32.63

*Three of the negative samples initially yielded undetermined Ct values for the IC and were reported as invalid. Repeat results were valid and negative.

Table: Clinical evaluation with nasopharyngeal samples.

SARS-CoV-2 concentration	Number of samples	Detection rate		Mean Ct		
		N	ORF1ab	N	ORF1ab	IC
2xLoD	20	20/20	20/20	38.01	37.77	31.98
4xLoD	20	20/20	20/20	37.12	36.32	32.11
10xLoD	1	1/1	1/1	35.46	34.72	31.64
20xLoD	1	1/1	1/1	35.46	34.23	32.13
50xLoD	1	1/1	1/1	33.27	32.92	29.86
100xLoD	1	1/1	1/1	31.78	31.43	30.46
200xLoD	1	1/1	1/1	32.95	31.49	32.08
250xLoD	1	1/1	1/1	31.85	30.49	32.04
500xLoD	1	1/1	1/1	30.40	29.73	30.24
Negative	94*	0/94	0/94	/	/	31.78

*One of the negative samples initially yielded undetermined Ct values for the IC and was reported as invalid. The repeat result was valid and negative.

As shown all positive samples at 2xLoD, 4xLoD, 10xLoD, 20xLoD, 50xLoD, 100xLoD, 200xLoD, 250xLoD and 500xLoD were positive and all negative samples were negative in the background of individual oropharyngeal swab and nasopharyngeal swab matrix.









References

1. *Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases: World Health Organization; 2020.*
2. *Innis, MA, et Al, and the PCR Protocols Applications: A Laboratory Manual , Academic, New York, 1989.*
3. *Mahony JB. Detection of Respiratory Viruses by Molecular Methods. Clinical Microbiology Reviews. 2008; 21 (4): 716-747.*
4. *China CDC Virus Disease Control and Prevention. Novel coronavirus nucleic acid detection primer and probe sequences (Specific Primers and Probes for Detection Novel coronavirus 2019) [EB / OL]., 2020-01-21.*

Revision history: Publication Number v1.0

Revision	Date	Description
1.0	March 20, 2020	New document

Key to symbols used

Symbol	Symbol Title and Reference Number
	Batch number
	Use-by date
	Temperature limit
	Contains sufficient for <n> tests
	Consult instructions for use
	Manufacturer
	This way up
	Fragile

For more information contact: COVID-19.TechnicalSupport@PerkinElmer.com



PerkinElmer, Inc.
7050 Burleson Road
Austin, Texas 78744
P: (800) 762-4000 or
(+1)203-925-4602
www.perkinelmer.com



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright ©2020, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc.
All other trademarks are the property of their respective owners.