

## New Tools in Food Analysis for Residue Screening and Authenticity

Jianru Stahl-Zeng, PhD Global Technical Manager – Market Development SCIEX



### **Today's Food Movement**

Emphasis on ensuring that food is safe and free of chemical residues and contaminants as it transits from the farm to the consumer.













Farming practices?

Growth climate and conditions? Adulteration?

#### Everything occurring between:

Manufacturing, processing, packaging, storage, transportation...



### **Target & Non-Target Screening on Food Testing**

### In 1 Injection.....





### So You Know What You Are Looking For?



But what could you be missing?



## LC-MS/MS with Multiple Reaction Monitoring (MRM)





### Triple Quad<sup>™</sup> 3500 System – Pesticide Screening

### Hundreds of Compounds in a Single Analysis



Over 400 MRM transitions monitored using the *Scheduled* MRM<sup>™</sup> Pro algorithm to quantify and identify pesticides in QuEChERS extract of fruit and vegetable samples at established maximum residue limits



### **QTRAP® 6500 System – Ion Ratios**

#### 0.1 ng/mL Avermectin B1a (10 µL injected)





## **Full Scan MS/MS using QTRAP® Technology**

Screening and Quantitation with Increased Confidence in Identification



### **Information Dependent Acquisition (IDA)**



### **Compound Identification using MS/MS Full Scan**



### Scheduled MRM-IDA-MS/MS complementary to MRM Ratios



Processing in MasterView<sup>™</sup> Software



### Why QTRAP<sup>®</sup> ? Increase Confidence in Identification

20 pesticides spiked into different fruit samples

Identification based on MRM ratio

- <section-header><section-header><text><text><list-item><list-item><list-item><list-item><list-item><list-item><list-item><list-item><list-item>
- 95% correct identification, 1% questionable (review), 4% false negative (outside acceptable MRM tolerance) – all confirmed by MS/MS library searching
- Identification based on MS/MS library searching
  - 99% correct identification, 1% questionable (review) confirmed with ratio



### **Pesticides in Yellow Pepper**



QuEChERS Phenomenex roQ (EN 15562 with PSA dSPE) Extract of 10 g with 10 mL and 10x dilution of extract with water



### **Pesticides in Spinach**



QuEChERS Phenomenex roQ (EN 15562 with PSA/GCB dSPE) Extract of 10 g with 10 mL and 10x dilution of extract with water

### Allergen Analysis by LC-MS/MS

- Over 150 million people worldwide suffer from a food allergy
  - Food allergy statistics in North America:
    - ~8% of children
    - 3.7% of adults
    - ... and rising
- What are allergens?
  - Complex mix of chemicals including proteins and sulfites that induce an adverse response in the body
- The cure?
  - There is none! The only solution is to avoid consumption of the allergen.
- As a result, accurate food labeling for allergens is extremely important and requires analysis.

![](_page_13_Picture_12.jpeg)

![](_page_14_Picture_1.jpeg)

### **Sample Preparation**

### Identification

![](_page_14_Figure_4.jpeg)

![](_page_14_Picture_5.jpeg)

### **Multi-Allergen Screening**

### Method built on the foundation of the work of Heick and Pöpping

![](_page_15_Picture_2.jpeg)

### Benefits of LC-MS/MS:

- Multiple allergens in one injection
- Fewer false positive or false negative results
- Multiple fragments for each allergen detected
- Full MS/MS
  fragmentation enables
  sequencing for highest
  level of selectivity
- Accurate quantitation

![](_page_15_Figure_9.jpeg)

Multi-allergen analysis

![](_page_15_Picture_11.jpeg)

#### More Confidence

### **Expending the Allergen-Method (1)**

- LC-MS/MS for nut allergen analysis:
  - Very large class
  - Includes tree nuts (i.e. pistachio, brazil, hazelnut, cashew, macadamia, walnut, pecan)
  - Includes peanut
  - Includes sesame

![](_page_16_Picture_6.jpeg)

![](_page_16_Picture_7.jpeg)

### Is what we eat authentic?

## Horse meat scandal highlights need for systematic testing

![](_page_17_Picture_2.jpeg)

By Joe Whitworth+ ■<sup>27</sup> 28-Aug-2013

Related tags: Brand protection, Eurofins, Danone, Horsemeat, , Testing, Regulation <u>Related topics</u>: Financial & Industry, Meat, fish and savoury ingredients

The horse meat scandal in Europe highlights the need for more systematic testing, according to Eurofins.

![](_page_17_Picture_6.jpeg)

Post a com

Post a com

## After horsemeat, pork may be next scandal, say farmers

By Fiona Barry 27, 27-Aug-2013

![](_page_17_Picture_9.jpeg)

### Authenticity of meat Horse, pork, duck, lamb, fox...

#### Food and Environme

#### an LC-MS/MS Be Used in Horse Meat Detection?

ephen Lock 3 SCIEX Warrington, Cheshire (UK)

#### verview

apid, robust, sensitive and specific LO-MG/MG assay has an developed for the simultaneous detection of horse meat at % levels in beef and the banned substance pheny/butazone TEI) using peptides markers for horse proteins and specific M transitions for BUTE.

#### troduction

wing the Food Standards Agency's (FSA) announcement in uary that horse and pig DNA had been identified in beef ducts sold by several supermarket chains, further testing oss Europe and beyond has revealed widespread incidence ch contamination.<sup>1</sup> However, most testing methods are ed on detection of species-specific DNA in meat, using the merase chain reaction (PCR) - which does not detect or of proteins. This is a concern because DNA can be easily upted or removed during standard meat processing and food ufacturing. As a result, horse tissue or other contaminants ain undetected in food samples, despite strong presence of contaminating proteins. An alternative protein-based method, 3A (enzyme-linked immunosorbent assay), can be used to plement DNA testing, but this method has limitations uding that it detects only one part of the protein and not iple protein markers.

E.C-MBUNG-based method presented offers a more accurate d reliable approach to meat speciation than PCR or ELIDAsed techniques or other indirect methods, and also allows for s detection of veterinary drug residues in the same analysis, ich is not possible by ELIDA or PCR.

e method was developed using an Exigent ekspert\* colo 200 UHCL system coupled with an AB SOLEX RAP<sup>6</sup> SSOL CLAMDMB system. The method uses multiple cition monitoring (IRMI) to detect peoplet markners for horse a la capabile of providing sequence information by scalining an hanced product ion (EPI) scan for each triggering MRM which he used to further confirm the peptide 1 / proteins and refore the species identity. This gives greater confidence for d testing when distinguishing between species; for example se and beef proteins may differ by as IBie as one or two ino acids.

![](_page_17_Picture_20.jpeg)

At the same time it is also possible to detect and quantity veterinary drug residues using the same extraction method an LC conditions by simply adding additional NRM transitions to to method. Here the nonsteroidal anti-inflammatory drug (NGAID) BUTE was detected in meats tamples.

#### Method Details

Standards

For the initial development work some of the target proteins we commercially available and therefore purchased as well as commercially available reference materials of port, beet, and horse meat and beet reference material which had been spiked at different levels with horse meat. A sample of lamb meat was obtained from a local supermarket.

A sigma standard of BUTE was not available at the time of thi work so BUTE had to be extracted from a sample of horse medicine.

Sample Preparation

The meat sample was homogenized using a food processor an mixed (2.g) with an extraction buffer containing this (2-amino-2hydroxymethyl-propane-1,3-diol), urea and acetonibrile (10 mL). The meat was broken up by shaking, uitra sonication (15 min) and agitated further using a rolicer mixer (45 min). This mixere

### Hot Topics and Trends – Top 10 Authenticity

![](_page_18_Picture_1.jpeg)

Related tags: Food fraud, Adulteration Related topics: Public Concerns, Testing

Dlive oil, fish and organic foods are at the highest risk of food fraud in Europe, accordi a new draft report from the European Union – but meat is not in the top ten, despite this rear's high-profile horse meat scandal.

![](_page_18_Picture_4.jpeg)

The report underlines that risk of food fraud is greatest when potential economic gains are large and chances of getting caught are slim. Most of the at-risk foods **RELATED NEWS:** 

Proposed regulation overhaul could cost industry, FSA warns Olive oil Meat and fish Organic produce Milk Grains Honey and maple syrup Coffee and tea Spices (chili, saffron) Fruit juices Wine

![](_page_18_Picture_9.jpeg)

### **Meat Speciation – Is it Authentic?**

- Similar workflow to allergen analysis can be used for meat speciation and authenticity testing
- Extraction, digestion, cleanup, LC-MS/MS

![](_page_19_Picture_3.jpeg)

![](_page_19_Picture_4.jpeg)

### **Sample Preparation**

![](_page_20_Figure_1.jpeg)

- Extraction of proteins
- Digestion into peptides
- Off-line SPE extraction
- LC-MS/MS analysis
  - Micro flow LC
  - Halo C18 50 x 0.5 mm
  - Gradient water/acetonitrile + 0.1% formic acid
  - 20 µL/min
  - SCIEX QTRAP<sup>®</sup> 5500 system
  - MRM-IDA-MS/MS (quantitation and identification)

![](_page_20_Picture_12.jpeg)

![](_page_20_Picture_13.jpeg)

### **Detection of Lamb, Beef, Horse Pork, and Horse Markers**

#### Meat Reference Materials Analyzed by LC-MS/MS

![](_page_21_Figure_2.jpeg)

Horse adulteration detectable at 1%

Updated method also used to identify duck, rat, fox... and to quantify veterinary drugs

![](_page_21_Picture_5.jpeg)

### **Detection of Pork Gelatin Markers**

![](_page_22_Figure_1.jpeg)

![](_page_22_Figure_2.jpeg)

Pork and beef gelatin markers detected in gummy bears, candy, and pharmaceutical capsule

![](_page_22_Picture_4.jpeg)

### Unknown Screening using High Resolution Accurate Mass LC-MS/MS

![](_page_23_Picture_1.jpeg)

![](_page_23_Picture_2.jpeg)

netalax	yl, subrange							meta	laxyl, subrange							
8e5	1		6.32						8e5 -							
6e5									6e5 -							
								da la								
465	1							-line	-H0 ]							
2e5									2e5 -		6.32					
0e0			<u> </u>						0:0							
	1 2 3	4 5 6	'	8 9 Time.mi	10 I		12 13 14 15 1	5 1/	1 2	3 4	5 6	/ 8 9 Time.n	i0 11	12 13	14 15	16
				1 800, 000												-
aster	View 📄 📩		X	<b>E</b> 1	4		<b>₩ \$</b>		? New Sessio	n	5	7				
	C T Wiff file Name	Sample Name	Number of positive	i	. *		1.111	Name	Formula	Intensity	Known Concentr.	Calculated Concentr.	Threshold (ratio of control)	Threshold (cps)	Control Intensity	la Di
1	Standard 2ppb	Sample 1	198	11	169	~	<b>VVVV</b>	metalaxyl	C15H21NO4	856968	10	119.0755	0.5	1000	71968	
~	Bok choy	Sample 1	1		151	~	<b>VVVV</b>	imidacloprid	C9H10CIN5O2	27009	10	19.7398	0.5	1000	13682	
5	Broccoli	Sample 1	6		248	$\checkmark$	<b>VVVV</b>	spirotetramat	C21H27NO5	180172	10	57.3822	0.5	1000	31399	
5	Cabbana	Samela 1	5		64	~	<b>VVVV</b>	cyprodinil	C14H15N3	11184	10	1.1189	0.5	1000	99960	
×.	Chinese berret	Comple 1			127	~	<b>VVV00</b>	fludicxonil	C12H6F2N2O2	2260	10	1.4983	0.5	1000	15086	
$\mathbf{x}$	Crimese broccou	Sample 1			243	$\checkmark$		spinetoram A	C42H69NO10	14065	10	75.2241	0.5	1000	1870	
~	Ciementine	Sample 1	2		245	~	<b>VV</b>	spinosyd A	C41H65NO10	1082	10	4.3323	0.5	1000	2497	
$\checkmark$	Kale	Sample 1	2		275	~	<b>VV</b>	triadimeton	C14H16CIN300	1024	10	0.6158	0.5	1000	16627	
$\checkmark$	Kohlrabi 2	Sample 1	1		143	~		furalaxyl	C17H19NO4	57449	10	122.6173	0.5	1000	4685	
$\checkmark$	Kohlrabi	Sample 1	3		250	~	<b>VAV00</b>	spiroxamine-1	C18H35NO2	502549	10	214.7741	0.5	1000	23399	
	Mustard green	Sample 1	0		110	~		etofenprox	C25H28O3	3154	10	1.6135	0.5	1000	19546	
~	Nappa cabbage	Sample 1	1		2	~	<b>VOV</b> 00	acephate	C4H10NO3PS	1324	10	5.3448	0.5	1000	2477	
	OJ 1	Sample 1	0		39	~		carboxin	C12H13NO2S	1647	10	0.2858	0.5	1000	57624	
	01.2	Sample 1	0		112	~		famoxadone	C22H18N2O4	1043	10	3.424	0.5	1000	3047	
	00000	Cample 1	2		249	~	<b>VOV0</b>	spiroxamine-1	C18H35NO2	502549	10	214.7741	0.5	1000	23399	
×	Unange	January 1			263	~	<b>VOV00</b>	tetramethrin-trans	C19H25NO4	26146	10	160.1716	0.5	1000	1632	
			50		262	~	<b>VOV00</b>	tetramethrin-cis	C19H25NO4	331.17836	+NH4	349.21218	0.01	349.21207	-0.1	_
ositive	result: equal or better															

![](_page_23_Picture_4.jpeg)

- Increased specificity through High Resolution acquisition for both MS and MS/MS data
  - Ability to screen for everything that can be ionised
- Calculate and postulate potential elemental formulas
  - Better to use both MS and MSMS information
- Accurate Mass MS/MS data, provides a chemical fingerprint of a molecule
  - Increases confidence in the data through less ambiguous results from MS only data (reduced number of false positive answers)
  - Allows the postulation of chemical structure based on fragmentation

![](_page_24_Picture_8.jpeg)

## **SCIEX TripleTOF® QTOF System**

![](_page_25_Picture_1.jpeg)

• Hi-Res accurate mass MS/MS analyser

© 2015 /

26

Sensitivity, speed and wide dynamic range

![](_page_25_Picture_4.jpeg)

### So, what is in my sample?

![](_page_26_Figure_1.jpeg)

#### **Targeted processing:**

#### Quickly screen against your target list:

- Identification using 4 criteria:
  - Retention time
  - Accurate MS
  - Isotopic patterns
  - Accurate MS/MS

#### Perform quantitation:

 Ability to quantitate how much is there if a compound is identified

#### Non-targeted processing:

#### Effectively find unknowns by:

- Sample-control comparisons
- Automatic formula finder
- Library searching
- ChemSpider search for structure elucidation

![](_page_26_Picture_17.jpeg)

![](_page_26_Picture_18.jpeg)

### Searching for what you know

### **Automatic Target Identification and Quantitation**

![](_page_27_Figure_2.jpeg)

### **Process Results Based on Targeted List with MS/MS**

![](_page_28_Figure_1.jpeg)

![](_page_28_Picture_2.jpeg)

### **Data EU Proficiency Test**

#### MS/MS Information is Crucial to Minimize False Positive Results

![](_page_29_Figure_2.jpeg)

High probability of false positive for Prometon without MS/MS information

![](_page_29_Picture_5.jpeg)

### **EU Proficiency Test Data - Results**

- Sample were analyzed using a SCIEX TripleTOF<sup>®</sup> 5600 system
- Samples were processed in MasterView<sup>™</sup> software
  - XIC table and library of 567 pesticides
  - Data processing < 30 seconds</li>
- 43 out of 44 pesticides in the provided leek and pear sample were correctly identified
  - One false negative result (low QuEChERS recovery for Ethoxyquin)
  - Mass errors between -1.7 and 2.1 ppm (positive polarity) and -0.7 and 2ppm (negative polarity) were well below required 5 ppm and retention time errors were well below required 2.5% (SANCO/12495/2011)
  - No false positive results due to MS/MS library searching

![](_page_30_Picture_9.jpeg)

### What's Really in my Sample? – Non Target Screening

![](_page_31_Figure_1.jpeg)

![](_page_31_Picture_2.jpeg)

### So, <u>WHAT</u> really is in my sample?

![](_page_32_Picture_1.jpeg)

Two streamlined data processing workflows:

![](_page_32_Picture_3.jpeg)

#### Non-targeted processing:

#### Effectively find unknowns by:

- Sample-control comparisons
- Automatic formula finder
- Library searching
- ChemSpider search for structure elucidation

![](_page_32_Picture_10.jpeg)

![](_page_32_Picture_11.jpeg)

## **Simplifying Data Complexity**

### **Comparative Unknown Screening** What is it and why is this a good approach?

Comparative screening = comparing your unknown sample to a known control sample to more quickly isolate unique features

vs. organic food product

vs. clean water sample

vs. un-treated sample

Comparative screening processing can reduce the list of compounds of interest from **tens of thousands** to **less than a hundred** 

> Get to the answer of what is in the sample faster and more accurately

![](_page_33_Picture_8.jpeg)

![](_page_33_Picture_9.jpeg)

#### TIC of Sample and Control

![](_page_34_Figure_2.jpeg)

![](_page_34_Picture_3.jpeg)

Intensity

### **Sample to Control Comparison**

![](_page_35_Figure_1.jpeg)

### **Automatic Formula Finding and MS/MS Search**

![](_page_36_Figure_1.jpeg)

1272 chromatographic peaks and 21 relevant signals

![](_page_36_Picture_3.jpeg)

### Non-Target Screening using MasterView™

### Review of Formula Finding and Library Searching Results

![](_page_37_Figure_2.jpeg)

![](_page_37_Picture_3.jpeg)

MS and MS/MS match plus ChemSpider hit count to quickly identify most likely formula

![](_page_37_Picture_5.jpeg)

MS

### **ChemSpider Search & MS/MS Interpretation**

![](_page_38_Figure_1.jpeg)

ChemSpider hits are automatically compared against MS/MS spectrum. SCIEX

### Hi-Res Accurate Mass MS/MS as a routine tool?

#### Making the 'Needle' in the 'Haystack' Easier to find

![](_page_39_Picture_2.jpeg)

![](_page_39_Picture_3.jpeg)

![](_page_40_Picture_0.jpeg)

Answers for Science. Knowledge for Life.™

![](_page_40_Picture_2.jpeg)

**PPCPs in Water** 

![](_page_41_Picture_0.jpeg)

### **Variable Window Calculator Tool**

Excel Sheet for Designing based on Constant Precursor Density

![](_page_42_Figure_2.jpeg)

- Open TIC of target Compounds (Pesticides) and extract out a MS spectrum of the whole LC run
- Generate m/z vs intensity list to paste into Excel
- Set assay parameters and build method
- www.absciex.com\VariableWindowsCalculator

Co o	omputation f Windows	AB SCIEX	
	SWATH Variable \	Nindow Assay Controls	
	Target number of windows:	(actual # may be less depending on min window setting)	
	Lower m/z limit: Upper m/z limit:	400 1500	
	Round bin edges to x figures:	1	
	Window overlap (Da)	1.0	
	Minimum window width (Da)	4	
	CES	5	

449.5	399.5
481.5	448.5
505.5	480.5
527.5	504.5
547.5	526.5
565.5	546.5
579.5	564.5
594.5	578.5
608.5	593.5
623.5	607.5
637.5	622.5
651.5	636.5
664.5	650.5
677.5	663.5
691.5	676.5
705.5	690.5
718.5	704.5
732.5	717.5
745.5	731.5
757.5	744.5
772.5	756.5

![](_page_42_Picture_9.jpeg)

### **Acquisition Method – Variable Mass Windows**

#### Example of Acquisition Method for a Mixture of 264 PPCP Compounds

![](_page_43_Picture_2.jpeg)

SCIEX

### **SWATH & IDA Acquisition**

#### Compounds can be identified for both SWATH and IDA data using MasterView

![](_page_44_Figure_2.jpeg)

![](_page_44_Picture_3.jpeg)

### **Results SWATH vs IDA**

#### Example Linuron, SWATH Window has more than one precusor

![](_page_45_Figure_2.jpeg)

![](_page_45_Picture_3.jpeg)

### **Results SWATH vs IDA**

#### Example Linuron, SWATH Window has more than one precusor

![](_page_46_Figure_2.jpeg)

![](_page_46_Picture_3.jpeg)

### **MS/MS SWATH vs IDA**

#### Intensity comparison

![](_page_47_Figure_2.jpeg)

![](_page_47_Picture_3.jpeg)

### **MS/MS SWATH vs IDA**

#### No MS/MS spectrum with IDA

![](_page_48_Figure_2.jpeg)

![](_page_48_Picture_3.jpeg)

- Non-targeted or unknown screening in food safety testing is possible
- Accurate Mass MS (but more importantly accurate mass MS/MS) data reduces matrix problems and increase selectivity for more confident identification
- Accurate mass MS/MS provides the capability of library searchable spectra that can be processed through ChemSpider for intuitive deconvolution
- Fast acquisition systems such as the TripleTOF allow all analytes ionised to be detected in both MS & MS/MS without reduction in data quality
- Software processing plays a big part....

![](_page_49_Picture_6.jpeg)

![](_page_50_Picture_0.jpeg)

# Answers for Science. Knowledge for Life.<sup>™</sup>

![](_page_50_Picture_2.jpeg)

- Andre Schrieber
- Jianru Stahl-Zeng
- Detlev Schleuder
- Harald Moeller
- Prof Amadeo Alba, EURL Almeria

![](_page_51_Picture_6.jpeg)