



Answers for Science.  
Knowledge for Life.™



## **Sciex X500R a new Screening Platform for Food, Environmental & Forensic Analysis**

Alexandre Paccou  
*Sr Manager, Support, EMEA*

For Research Use Only. Not for use in diagnostic procedures.

RUO-MKT-11-3046-A

# What's New?

## SCIEX X500R QTOF System powered by SCIEX OS Software



SCIEX X500R QTOF system

SCIEX ExionLC™ AC system

SCIEX OS software

# Outline

- Hardware
  - SCIEX ExionLC™ Systems
    - Fully controlled by SCIEX OS software
    - Improved software integration for better stability
  - SCIEX X500R QTOF System
    - N-optic design
    - Heated TOF path
    - Minimized footprint, engineered for simplicity and service accessibility
- Software
  - SCIEX OS Software
    - New user interface
    - Simultaneous identification and quantitation
    - Automatic unknown identification
- Application data
  - Target identification and quantitation
  - Unknown screening



# Introducing the SCIEX X500R QTOF System

## Design Improvements and Details

### Engineered for simplicity

Optimal performance delivered through adjustment of only 2 voltage variables.

### TwinSpray

An independent calibrant delivery path for reliable auto-calibration.

### Service accessibility


Easy QJet® and Turbo pump access for fast and efficient maintenance, increasing system uptime.

### Minimized footprint

The benchtop stature (110 x 57 x 112 cm)\* occupies less lab space than any other HRMS system on the market.


### Integrated calibration

Maintains mass calibration through long runs without effect on sample flow.



# Introducing the SCIEX X500R QTOF System

## Design Improvements and Details





# Introducing the SCIEX X500R QTOF System

## Design Improvements and Details

### Dynamic Background Calibration

Use of background ions to maintain stable mass axis without calibrant.

### TwinSprayer






### TwinSpray Probe



# History 5

## TOF Details

- Simplified ion transfer optics
- No slits... all ions go in (increased sensitivity)
- Focusing... high resolution without compromise in sensitivity
- Single pulse acceleration (6kV)
- 2 mirrors and direct thermal stabilization



# Perfect Balance to Elevate Your Lab's Performance

## SCIEX X500R QTOF System




- The first robust, high performance high resolution MS/MS system designed for routine use.
  - Sensitivity to easily detect required drug cutoff concentration levels
  - Resolving power to remove interference from complex forensic matrices
  - Linearity to quantify up to 4 orders of magnitude
  - Identify compounds based on mass accuracy
  - Confident identification based MS/MS (IDA, MS/MS<sup>ALL</sup> with SWATH<sup>TM</sup> Acquisition and ion ratios)
  - Industry leading robustness of Turbo V<sup>TM</sup> source and Curtain Gas<sup>TM</sup> interface




# The SCIEX X500R QTOF System

Sensitivity and Resolution (20 ng/mL forensic drugs in urine, 10  $\mu$ L injected)



# Resolution


Resolution from Dec 9 to Dec 16



TOF-MS

TOF-MS/MS

# The SCIEX X500R QTOF System – Linearity




Linearity (Asenapine 0.5 to 1000 ng/mL)

RUO-MKT-11-3046-A



# Mass Accuracy in urine samples without autocalibration


ppm in all injections over a 10-hour period (no calibration in between)



**This is a plot for 6 urine samples with same analytes in different concentrations over a 10-hour period without any calibration.**

# MS/MS<sup>ALL</sup> using SWATH<sup>®</sup> Acquisition

A Mode of Data Independent Acquisition Providing all MS/MS Spectra  
Wide Q1 Isolation Window is Stepped Across the Mass Range




SWATH<sup>®</sup> acquisition can be used with variable Q1 isolation windows across the mass range




## SWATH<sup>®</sup> Acquisition

- A mode of Data Independent Acquisition
- Provides MS/MS<sup>ALL</sup>

Regular MS/MS




SWATH MS/MS



# The SCIEX X500R QTOF System – SWATH® Acquisition


20 - "3,4-dimethylmethcathinone 1" (Standard) 192.1283 - 19...ning\Data\DOA SWATH290715.wiff2, (sample Index: 8)  
Area: 1.908e5, Height: 2.781e4, RT: 4.41 min



▼ Peak Details

Precursor m/z	Retention Time (min)	Ion Ratio
192.138	4.41	1.0000


● Spectrum from DOA SWATH290715.wiff2 (sample 1...F MSMS [189 - 200] (30 - 500) from 4.415 min  
subtracted by: [Spectrum from DOA SWATH2907...[189 - 200] (30 - 500) from 4.08 to 4.100 min]  
subtracted by: [Spectrum from DOA SWATH2907...[189 - 200] (30 - 500) from 4.14 to 4.71 min]



▼ Library Search Results

Name	CAS#	Formula	MW (Da)	Fit	Rev. Fit	Purity	CE (eV)
3,4-Dimethylmethcathinone	1082110-00-6	$C_{12}H_{17}NO$	191.27327	92.9	100.0	92.9	35


20 - 4-ethylmethcathinone 1 (Standard) 192.1283 - 19...ning\Data\DOA SWATH290715.wiff2, (sample Index: 8)  
Area: 2.031e5, Height: 2.429e4, RT: 4.61 min



▼ Peak Details

Precursor m/z	Retention Time (min)	Ion Ratio
192.138	4.61	1.0000

● Spectrum from DOA SWATH290715.wiff2 (sample 1...F MSMS [189 - 200] (30 - 500) from 4.612 min  
subtracted by: [Spectrum from DOA SWATH2907...[189 - 200] (30 - 500) from 4.234 to 4.514 min]  
subtracted by: [Spectrum from DOA SWATH2907...[189 - 200] (30 - 500) from 4.793 to 5.073 min]



▼ Library Search Results

Name	CAS#	Formula	MW (Da)	Fit	Rev. Fit	Purity	CE (eV)
4-Ethylmethcathinone	1225622-14-9	$C_{12}H_{17}NO$	191.27327	98.3	100.0	98.3	35

Confident Identification of Isomers using MS/MS

# Acquisition method

## SWATH with variable windows, 8 min

### TOF MSMS

TOF start mass  Da TOF stop mass  Da Dynamic collision energy ☐

Accumulation time  s Charge state

### Mass Table

[Autofill SWATH windows...](#)

	Precursor ion start mass (Da)	Precursor ion stop mass (Da)	Declustering potential (V)	DP spread (V)	Collision energy (V)	CE spread (V)
1	120.0000	140.0000	60	0	35	15
2	139.0000	165.0000	60	0	35	15
3	164.0000	195.0000	60	0	35	15
4	194.0000	205.0000	60	0	35	15
5	204.0000	218.0000	60	0	35	15
6	217.0000	240.0000	60	0	35	15
7	239.0000	258.0000	60	0	35	15
8	257.0000	273.0000	60	0	35	15
9	272.0000	290.0000	60	0	35	15
10	289.0000	304.0000	60	0	35	15
11	303.0000	313.0000	60	0	35	15
12	312.0000	326.0000	60	0	35	15
13	325.0000	333.0000	60	0	35	15
14	332.0000	358.0000	60	0	35	15
15	357.0000	450.0000	60	0	35	15
16	449.0000	500.0000	60	0	35	15

# Acquisition method

## SWATH with variable windows, 8 min

### TOF MSMS

TOF start mass  Da

TOF stop mass  Da

Dynamic collision energy ☐

Accumulation time  s

Charge state

### Mass Table


[Autofill SWATH windows...](#)

	Precursor ion start mass (Da)	Precursor ion stop mass (Da)	Decustering potential (V)	DP spread (V)	Collision energy (V)	CE spread (V)
1	120.0000	140.0000	Amphetamine (136)			
2	139.0000	165.0000	Methamphetamine (150)			
3	164.0000	195.0000	MDA (180)			
4	194.0000	205.0000	MDMA (194)	Nortapentadol (204)		
5	204.0000	218.0000		Tapentadol (222)		
6	217.0000	240.0000	Norketamine (224)			
7	239.0000	258.0000	Ketamine (238)	Nortramadol (250)		
8	257.0000	273.0000		Tramadol (264)	Norvenlafaxine (264)	Nordoxepin (266)
9	272.0000	290.0000	Morphine (286)	Norchlordiazepoxide (286)	Verlafaxine (278)	Doxepin (280)
10	289.0000	304.0000	Codeine (300)	Chlordiazepoxide (300)	Norsertaline (292)	Norfluoxetine (296)
11	303.0000	313.0000	Alprazolam (309)	Norcitalopram (311)	Sertraline (306)	Fuoxetine (310)
12	312.0000	326.0000	Hydroxyalprazolam (325)	Citalopram (325)	Norclozapine (313)	THC (315)
13	325.0000	333.0000	Norpropoxyphene (326)		Clozapine (327)	THC-OH (331)
14	332.0000	358.0000	Propoxyphene (340)	Triazolam (343)		
15	357.0000	450.0000	Norbuprenorphine (414)	Hydrotriazolam (359)	Norverapamil (441)	
16	449.0000	500.0000	Buprenorphine (468)		Verapamil (455)	



# The SCIEX X500R QTOF System – MRM<sup>HR</sup>

Increased selectivity with MRM<sup>HR</sup>: Avoiding false negatives  
Feed Sample Tested Positive for NP-Semicarbazide





# Acquisition method

*Scheduled MRM<sup>HR</sup>, 8 min, Fragment  $\pm 10$  m/z*

## TOF MS

TOF start mass	100 Da	Declustering potential	60 V	Collision energy	10 V
TOF stop mass	1000 Da	DP spread	0 V	CE spread	0 V
Accumulation time	0.1 sec				

## TOF MSMS

Mass Table


☒ Apply fragment ion mass
 ☐ Apply TOF start/stop mass
 ☒ Apply Scan Schedule
 [Import and autofill...](#)
[Sort by precursor ion](#)

	Compound ID	Precursor i...	Fragment ion (Da)	Accumul...	Declus...	Collision...	Retention ti...	Retention time tolerance (+/- sec)
1	6-MAM	328.15	211.0747	0.0100	60	33	2.04	15
2	7-Aminoclonazepam	286.07	222.1025	0.0100	60	32	2.67	15
3	7-Hydroxymisragline	415.22	190.0864	0.0100	60	40	2.96	15
4	Acetyl Fentanyl	323.21	105.0700	0.0100	60	48	3.00	15
5	Alpha-Hydroxyalprazolam	325.09	297.0665	0.0100	60	31	4.18	15
6	Alpha-hydroxymidazolam	342.08	203.0377	0.0100	60	34	4.51	15
7	Alpha-hydroxytriazolam	359.05	331.0272	0.0100	60	35	3.95	15
8	Alpha-PPP	204.14	105.0699	0.0100	60	35	2.04	15
9	Alpha-PVP	232.17	161.0963	0.0100	60	20	2.61	15
10	Alprazolam	309.09	281.0730	0.0100	60	35	4.49	15
11	AM-2201 4-OH pentyl	376.17	155.0492	0.0100	60	35	6.07	15
12	Amitriptyline	278.19	117.0702	0.0100	60	31	4.53	15
■ ■ ■								
		Precursor	+/- 10 MS/MS	100 Hz			RT	RT half window
84	Zolpidem	308.18	235.1262	0.0100	60	45	3.06	15
85	THC-COOH	345.21	299.2157	0.0100	60	30	6.34	15


# The SCIEX X500R QTOF System – MRM<sup>HR</sup>

*Scheduled MRM<sup>HR</sup> (vs TOF scan)*

Buphedrone (5 ng/mL in urine, 10-fold dilution, 10  $\mu$ L injection)




**MRM<sup>HR</sup>: 178.1  $\rightarrow$  131.070  $\pm$  0.01  $m/z$**



**TOF-MS: 178.1226  $\pm$  0.005  $m/z$**

# The SCIEX X500R QTOF System – MRM<sup>HR</sup>


Scheduled MRM<sup>HR</sup>, 372.2→169.0644±0.0100 m/z



JWH-122 5-OH pentyl in urine (Urine was diluted 10-fold, 10  $\mu$ L injection)

# The SCIEX X500R QTOF System – Performance


Sensitivity and Resolution (1 ng/mL Pesticides , 5  $\mu$ L injected, XIC  $\pm$  5 mDa)





# The SCIEX X500R QTOF System – Performance

## Linearity and Mass Accuracy (Paclobutrazol 0.1 to 10,000 ng/mL)




Mass error = -0.2 to 0.91 ppm, resolution of 30500 @ 294



# The SCIEX X500R QTOF System – Performance


## MS/MS Identification of Thiachloprid (10 ng/mL)



Review of XIC, MS and MS/MS

# Information Dependent Acquisition of MS/MS (IDA)


Provides MS/MS Spectra with High Selectivity (Q1 Resolution unit)



IDA can be used with and without inclusion list for target or non-target screening

# The SCIEX X500R QTOF System – IDA


## Confident Identification of Isomers using MS/MS (Prometon vs. Terbumeton)



Identification based on MS/MS library searching (Q1 at unit resolution)

# MRM<sup>HR</sup>

## Combining Selectivity of MRM and High Resolution




Unit Q1 isolation, fragmentation using optimized CE and narrow XIC of product ion



# The SCIEX X500R QTOF System – MRM<sup>HR</sup>

## 24 µg/L Glyphosate in a Beer Sample






# Perfect Balance to Elevate Your Lab's Performance

## SCIEX X500R QTOF System Powered by SCIEX OS Software




- The first robust, high performance high resolution MS/MS system designed for routine use.
  - Sensitivity to easily detect required drug cutoff concentration levels
  - Resolving power to remove interference from complex forensic matrices
  - Linearity to quantify up to 4 orders of magnitude
  - Identify compounds based on mass accuracy
  - Confident identification based MS/MS (IDA, MS/MS<sup>ALL</sup> with SWATH™ Acquisition and ion ratios
  - Industry leading robustness of Turbo V™ source and Curtain Gas™ interface
- Simple software workflows that deliver reliable results, it's the solution that's ready to meet your challenges today and for the future, from the trusted LC-MS/MS industry leader.

# Targeted Data Processing Workflow




XIC generation for each target



Identification



Quantitation



RT ✓

Mass

✓


MS/MS search ✓

# Identification and Quantification in SCIEX OS Software

## 3. Review Quantitative and Qualitative Results

Traffic Lights to indicate confidence

Flagging when above cutoff concentration




Review of peak integration, spectra and calibration line



# Identification and Quantification in SCIEX OS Software

## 3. Review Quantitative and Qualitative Results




Targeted Data Processing of TOF-MS -Information Dependent Acquisition- TOF-MS/MS acquired data



# Identification and Quantification in SCIEX OS Software

Highest confidence gained through library MS/MS comparison




XIC chromatograms

MS spectra

MS/MS spectra

# Confirmation and quantitation using SWATH® acquisition for forensic drug screening



1. Identification based on unique fragment ions and their ratios as well as MS/MS library searching
2. Quantitation on both precursor ion and fragment ions

# SCIEX OS Software – Unknown Identification




## Software Improvements and Details

**Automatic sample-control comparison**

**Non-Target peak finding**

**ChemSpider searching**


Sorting by number of references and automatic structure elucidation to identify unknowns




**Automatic MS/MS library searching and formula finding**



# Non-Targeted Data Processing Workflow




Sample-control comparison



Mass and RT

Identification




Formula finding ✓

MS/MS search ✓



# Unknown Identification in SCIEX OS Software


## 1. Automatic Sample-Control Comparison and Compound Identification



Automatic formula finding

# Unknown Identification in SCIEX OS Software

## 1. Automatic Sample-Control Comparison and Compound Identification



# Summary

- Hardware
  - SCIEX ExionLC™ Systems
    - Fully controlled by SCIEX OS software
    - Improved software integration for better stability
  - SCIEX X500R QTOF System
    - N-optic design
    - Heated TOF path
    - Minimized footprint, engineered for simplicity and service accessibility
- Software
  - SCIEX OS Software
    - New user interface
    - Simultaneous identification and quantitation
    - Automatic unknown identification
- Application data
  - Target identification and quantitation
  - Unknown screening



Thank you for your attention!



You ready to get behind the wheel?



# Trademarks/Licensing

For Research Use Only. Not for use in diagnostic procedures.

The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license. AB Sciex is operating as SCIEX.

© 2016 AB Sciex.