Avio ICP-OES applications



Helmut Ernstberger 13.6.2019



Wastewater Analysis with the Avio® 500 ICP-OES and EPA Method 200.7





Introduction

- Pollution control & prevention is critical for human and environmental health
- Water must be monitored at all stages
- U.S EPA Method 200.7
 - Determination of metals & trace elements in waters and wastes by ICP-OES
 - Variety of implementations
- Focus on wastewaters
 - High elemental concentrations & difficult matrices \rightarrow perfect for ICP-OES

Sample & Sample Preparation

- Samples: Wastewater Reference Materials
 - Wastewaters C & D (High Purity Standards)
 - Wastewaters Low & High (Inorganic Ventures)
- Sample Preparation
 - Prepare as directed on reference materials
 - No further dilution
- Calibrations

0

External	Element	Standard 1 (mg/L)	Standard 2 (mg/L)	Standard 3 (mg/L)
	Ag, Al, As, B, Ba, Be, Cd, Ce, Co, Cr, Cu, Fe, Li, Mn, Mo, Ni, P, Pb, Sb, Se, Si, Sn, Sr, Ti, Tl, V, Zn	0.5	1	5
	Na, Mg, K, Ca	20.5	51	105

Internal Standards: Sc & Y



Instrumental Conditions: Avio® 500 ICP-OES

Instrumental Parameters

Parameter	Value		
Nebulizer	MEINHARD® Type K		
Spray Chamber	Baffled glass cyclonic	1	
Sample Uptake Rate	1.0 mL/min		
RF Power	1500 W		
Injector	2.0 mm id Alumina		
Nebulizer Gas Flow	0.70 L/min		
Auxiliary Gas Flow	0.2 L/min		
Plasma Gas Flow	8 L/min	110113	
Integration Range	1-10 sec		
Sample Uptake Tubing	Black/Black (0.76 mm id)		
Internal Standard Tubing	Green/Orange (0.38 mm id)		
Drain Tubing	Red/Red (1.14 mm id)		
Replicates	4		



Instrumental Conditions: Avio® 500 ICP-OES



Elements, Wavelengths, and Plasma View

Element	Wavelength (nm)	Plasma View	Element	Wavelength (nm)	Plasma View
Ag	328.068	Radial	Мо	203.845	Axial
Al	394.401	Radial	Na	589.592	Radial
As	188.979	Axial	Ni	231.604	Axial
В	249.677	Axial	Р	178.221	Axial
Ва	493.408	Radial	Pb	220.353	Axial
Be	313.107	Axial	Sb	206.836	Axial
Ca	315.887	Radial	Se	196.026	Axial
Cd	214.440	Axial	Si	251.611	Axial
Со	228.616	Axial	Sn	189.927	Axial
Cr	267.716	Axial	Sr	421.552	Radial
Cu	324.752	Axial	Ti	334.940	Axial
Fe	238.204	Radial	TI	190.801	Axial
K	766.490	Radial	V	292.402	Axial
Li	670.784	Radial	Zn	206.200	Axial
Mg	285.213	Radial	Sc (int std)	361.383	Radial
Mn	257.610	Axial	Y (int std)	371.029	Axial

Results: Linear Dynamic Range

 Highest concentration that will recover within 10% based when analyzed against a normal calibration curve

Elements	Linear Range (mg/L)
Al, Be, Mn, Ti	50
Ag, As, B, Ba, Cd, Co, Cr, Cu, Fe, Li Mo, Ni, P, Pb, Sb, Se, Si, Sn, Sr, Tl, V, Zn	100*
Mg	200
Ca, K, Na	500*

- * = highest concentration tested
- Options for extending the linear range
 - Use a less sensitive sample introduction system
 - Vary the torch position

- Select less sensitive wavelengths
- Change the viewing height in the plasma



Results: Initial Quality Control

- Initial Performance Check (IPC)
 - Made from same stock solutions as the calibration standards
 - Analyzed after the calibration standards
 - All elements must recover within 5% of their true concentration
- Quality Control Sample (QCS)
 - Made from 2nd source standards
 - Same concentration as the IPC
 - Analyzed after the IPC
 - All elements must recover within 5% of their true concentration



Results: Initial Quality Control

All Elements Recover within 5%



Results: Method Detection Limits

 3.14 * Standard Deviation of 7 measurements of a blank spiked at 2-3 times the instrument detection limit



MDLs much lower than low reference material



Results: Spectral Interference Check

- Run 300 mg/L Fe + 200 mg/L AI
 - Specified in 200.7
- Look for false positives on analytes
- Confirm that false positives result from contamination and not interferences
- Spectral resolution of Avio 500 allows elements to be measured without the need for IECs



Results: Accuracy



• Run wastewater reference materials

Recoveries within 10% for all elements in 4 different reference materials



Results: Accuracy

- Elements in EPA 200.7 not certified in the reference materials
 - B, Li, P, Si \rightarrow Spike at 0.6 and 1.2 ppm
 - Na, Mg, Ca, K \rightarrow Spike at 75 and 150 ppm



Recoveries within 10% for all elements in 4 different RMs



Results: Stability

- Run IPC every 10 samples during a 10-hour wastewater analysis
 - Recoveries must be within 10% of true values



All recoveries within 10% over 10 hours



Results: Stability

- Monitor internal standards to check for instrument drift
 - Sc in radial mode, Y in axial mode



Variations of 5% over 10 hours



Summary

- The Avio[®] 500 ICP-OES easily analyzes wastewaters following EPA Method 200.7
- Take advantage of both axial and radial plasma viewing
- Accurate and stable analyses
- Vertical torch and low argon consumption provide a robust, cost-effective way to analyze wastewaters



To download the full application note, visit www.perkinelmer.com/avio 500



Analysis of In-Service Oils Following ASTM D5185 with the Avio 200 ICP-OES





Introduction

- With heavy machinery, important to monitor its status to prevent breakdown and costly repairs
- Monitor oil & lubricants for metal content
 - Engines, transmissions, gearboxes, and other areas
- Metal content provides information about the conditions of the engine and oil itself
 - Diagnostic for when oil needs changing
- Important elemental concentrations > 1 ppm
 - ICP-OES is ideal technique



Introduction: ASTM D5185 for In-Service Oils

Parameter	Specified by D5185	Common Implementation
Elements	Al, Ba, B, Ca,Cr, Cu, Fe, Pb, Mg, Mn, Mo, Ni, P, K, Na, Si, Ag, S, Sn, Ti, V, Zn	Elements important to the components being tested
Sample Preparation	By weight	By volume
Internal Standard	Cd, Co, or Y	Со
QC Frequency	Every 5 samples	Varies between 10-25 samples
QC Limits	± 5%	± 10%

- For high-throughput labs, simultaneous ICP-OES preferred
- For low-throughput labs, the Avio 200 is ideal

Samples & Sample Preparation

- Samples
 - In-Service Oils
- Sample Preparation
 - Dilute10x (by volume) with V-Solv spiked with 40 ppm Co (internal standard)
- Calibrations
 - Blank = 75 cSt base oil
 - Calibration Standards
 - V-23: 50, 100, 500 ppm
 - MA4 for Ca, P, Mg, Zn
- QC
 - 50 ppm V-23 + MA4



Instrumental Conditions: Avio[™] 200 ICP-OES

Elements and Wavelengths

Element	Wavelength (nm)		
Al	394.401		
Ag	328.068		
В	249.677		
Ba	232.527		
Ca	315.887		
Cd	228.802		
Cr	205.560		
Cu	324.752		
Fe	259.939		
К	766.490		
Mg	279.077		
Мо	203.845		
Mn	257.610		
Na	588.995		
Ni	232.003		
Р	214.914		
Pb	220.353		
Sb	217.582		
Si	288.158		
Sn	189.927		
Ti	334.940		
V	292.464		
Zn	213.857		
Co (internal standard)	228.616		

Instrumental	Parameters
--------------	------------

Parameter	Value	
Nebulizer	GemCone™ (modified Babbington)	
Spray Chamber	Baffled glass cyclonic	
RF Power	1500 W	
Torch	3-Slot Avio torch for organics	
Injector	1.2 mm ceramic	
Plasma Gas Flow	10 L/min	
Aux Gas Flow	0.8 L/min	W Argon
Nebulizer Gas Flow	0.40 L/min	FIOWS
Torch Position	-4	
Sample Uptake Rate	2.0 mL/min	
Sample Uptake Tubing	Black/Black (0.76 mm id), Viton	
Drain Tubing	Red/Red (1.14 mm id), SolvaFlex	
Read Delay	18 sec Short In	tegration
Replicates	2 Tin	nes
Rinse Between Samples	12 sec (fast pump at 6 mL/min)	
Auto Integration Range	0.2 - 1.0 sec	

Instrumental Conditions: Avio[™] 200 ICP-OES





- Correct Plasma Position
 - Carbon "bullet" just below 2nd flat plate



Results: Typical In-Service Oil Sample

Element	Concentration (ppm)	RSD (%)
Al	16	3.1
Ag	< 0.5*	-
В	81	2.8
Ba	< 0.5*	5.9
Ca	2441	0.21
Cd	< 0.5*	-
Cr	1	2.6
Cu	1	0.26
Fe	27	0.44
K	5	8.0
Mg	17	2.2
Mo	25	0.87
Mn	1	0.69
Na	236	0.44
Ni	2	12
Р	847	1.9
Pb	7	0.65
Sb	< 0.5*	-
Si	9	4.2
Sn	< 0.5*	-
Ti	< 0.5*	-
V	< 0.5*	-
Zn	974	2.8

* = lower quantification limit



Results: Final QC Sample after 20 Sample Analytical Run

	Final QC		
Element	Concentration (ppm)	% Recovery	
Al	50	100	
Ag	48	96	
В	54	108	
Ba	49	98	
Ca	45	90	
Cd	50	100	
Cr	50	100	
Cu	50	100	
Fe	53	106	
К	49	98	
Mg	46	92	
Mo	49	98	
Mn	49	98	
Na	49	98	
Ni	49	98	
Р	46	92	
Pb	50	100	
Sb	50	100	
Si	48	96	
Sn	49	98	
Ti	50	100	
V	49	98	
Zn	47	94	

• All recoveries within 10%

Summary

- The Avio 200 ICP-OES easily analyzes in-service oils following a common implementation of ASTM method D5185
 - Appropriate for low-volume labs
- Benefits of using the Avio 200 ICP-OES
 - Using only 8 L/min of argon, significant cost savings are realized
 - Vertical torch allows excess oil to drain
 - Short integration times allow for faster analyses



To download the full application note, visit www.perkinelmer.com/avio200



Analysis of SiO₂- and TiO₂-Containing Medications for Class 1 & 2A Elements with ICP-OES Following USP 232/233 Guidelines





Introduction: <USP 232> Defines Element & Impurity Limits





Introduction:

<USP 233> Defines Sample Prep & Analysis Requirements

Sample Prep

- Closed-vessel digestion
- Analysis
 - Accuracy
 - Spikes at 50%, 100% & 150% of maximum PDEs
 - Recoveries must be 70%-150%

• Repeatability

- Measure 6 independent samples spiked at 100% of target limit
- % RSDs over 6 samples must be < 20%

Ruggedness

- Measure 6 solutions on different days, with different instrument, or with different analyst
- %RSDs over 12 measurements must be < 25%

System Suitability

- Measure high standard at beginning & end of each batch
- Difference between results must be < 20%

Introduction: Samples

- Pharmaceutical tablets may contain a wide variety of excipients
 - Most are easy to digests
 - SiO₂ and TiO₂ are common and difficult to digest
- SiO₂ and TiO₂ are common excipients in a variety of over-thecounter medications
 - Acid reducers
 - Sleep aids
 - Motion sickness medication



Goal

- Develop sample preparation & analytical methodology for the analysis of Class 1 and 2A elements in SiO₂- and TiO₂containing tablets with the Avio 200
- Class 1 & 2A elements
 - As, Cd, Co, Hg, Ni, Pb, V
 - Most common to analyze
- SiO₂ & TiO₂ medications
 - Challenging sample preparation
- Avio 200
 - Ideal for smaller number of elements
 - Avio 500 could also be used

Samples

All samples purchased from local stores

Medication Type	Excipient	Active Ingredient	Daily Dose	Mass Per Tablet (g)
	TiO ₂	Ranitidine, 75 mg	2 tablets	0.13
Acid Reducer		Ranitidine, 150 mg	2 tablets	0.32
		Famotidine, 20 mg	2 tablets	0.21
Cloop Aid	SiO ₂	Diphenylhydramine, 25 mg	2 tablets	0.43
Sleep Alu		Diphenylhydramine, 25 mg	2 tablets	0.42
Motion Sickness	SiO2	Meclizine, 25 mg	2 tablets	0.20

- 3 brands of acid reducers
- 2 brands of sleep aids
- 1 brand of motion sickness



Sample Preparation: Closed Vessel Microwave Digestion

Titan MPS

Sample Amounts & Acids per Digestion Vessel

Sample	HNO ₃	HCl	HF	Water
	70% (mL)	35% (mL)	49% (mL)	(mL)
Tablet	1.5	0.5	0.5	7.5

• HNO₃

- Keeps elements in solution
- HCI
 - Specifically for Hg
- HF
 - $^\circ~$ Dissolve SiO_2 and TiO_2
- Water
 - Keep acid concentrations to a minimum
 - 3% HNO₃ + 1% HCl in final solutions



Sample Preparation: Closed Vessel Microwave Digestion





Titan Digestion Program

Step	Temperature (°C)	Pressure (Bar)	Ramp (Min)	Hold (Min)	Power (%)
1	160	35	5	1	90
2	190	35	5	5	100
3*	50	35	1	15	0

- Digestion Time = 16 minutes
- Cooling = 16 minutes



Sample Preparation: Closed Vessel Microwave Digestion

Complex HF

Titan MPS

- Allows use of standard sample introduction system
 - Can be skipped if use HF-resistant sample intro system
- Add 3 mL of saturated boric acid solution to each vessel
- Digest with following program

Titan Digestion Program

Step	Temperature (°C)	Pressure (Bar)	Ramp (Min)	Hold (Min)	Power (%)
1	190	35	5	10	90
2*	50	35	1	15	0



Calibration

- Calibration standards based on the J-Value
 - Defined in USP <233>
- Calibrate with 0.5J and 1.5 J standards
 - External calibration curves
 - Specified in USP <233>

Calibration standards made in 3% HNO₃ + 1% HCI

Element	0.5J (mg/L)	1J (mg/L)	1.5J (mg/L)
Cd	0.025	0.05	0.075
Pb	0.025	0.05	0.075
As	0.075	0.15	0.225
Hg	0.15	0.3	0.45
Со	0.25	0.5	0.75
V	0.50	1	1.5
Ni	1	2	3

Instrumental Conditions: Avio 200

Parameter	Value
Nebulizer	MEINHARD® Type K, glass
Spray Chamber	Baffled glass cyclonic
Sample Uptake Rate	1.0 mL/min
RF Power	1500 W
Injector	2.0 mm id Alumina
Nebulizer Gas Flow	0.70 mL/min
Auxiliary Gas Flow	0.2 L/min
Plasma Gas Flow	8 L/min
Torch Position	-3
Plasma View Mode	Axial
Replicates	3

Standard Sample Intro Components & Conditions

Element	Wavelength (nm)
As	193.696
Cd	214.440
Со	238.892
Hg	194.168
Ni	231.604
Pb	220.353
V	309.310
Y (Int std)	371.029

Results: Evaluation of Interferences

- 1.5 J single element standards of all analytes analyzed
- Overlay spectra to look for interferences on each analyte
- No interferences observed



Results: System Suitability

- USP <233> System Suitability Criteria (i.e. Drift Check)
 - Difference between 1.5J standard run at the beginning & end of a batch must be < 20%



System Suitability: Pass



Results: Accuracy

- USP <233> Accuracy Criteria
 - Pre-digestion spike recoveries at 0.5J, J, &1.5J must be between 70-150%



Recoveries within 10% for all spikes of all analytes in TiO₂ and SiO₂ Medications

Accuracy: Pass



Results: Repeatability

- USP <233> Repeatability Criteria
 - Same medication must be prepared 6 times with J pre-digestion spikes
 - %RSD of the measurements for these 6 samples must be < 20%



RSDs < 5% for all spikes of all analytes in TiO_2 and SiO_2 Medications

Repeatability: Pass



Results: Ruggedness

- USP <233> Ruggedness Criteria
 - Measure 6 samples used in Repeatability study on 2 different days
 - %RSD of the 12 measurements must be < 25%



RSDs < 5% for all spikes of all analytes in TiO_2 and SiO_2 Medications

Ruggedness: Pass



Results: Sample Analysis

- Concentrations of all elements in all samples < 0.3J
- 0.3J is considered an actionable level

Summary

- Avio 200 meets USP <233> criteria for Class 1 and 2A elements in medications containing TiO₂ or SiO₂ as excipients
 - Class 1 & 2A elements are the most commonly analyzed – perfect for Avio 200
 - In low throughput labs, the Avio 200 can also be used for Class 2B and 3 analytes
- Closed-vessel microwave digestion provides complete digestion with minimal acid
- Results easily surpass USP <233> criteria for accuracy, repeatability, ruggedness, and system suitability



Analysis of SiO₂- and TiO₂-Containing Medications Using ICP-OES Following USP 232/233 Guidelines with Software Designed to Aid in 21 CFR Part 11 Compliance

introduction

Since early 2016, manufactures of phermanucical products are mandated to comply with USP =2012-2210 - requirements around the analysis of elemental inguithment in medications.¹⁴ A detailed description of the explorement of USP =2222-2230

Is available, for only a brief description is given here. Betweets are categorized into four decises, with Class 1 long the most track and Class 20 elements being likely contravis and how the nanoullabeling particulations (i.e. from mixing exploreer, stables stead vessel), etc.). Therefore, these are the most important devisers to measure in pharmanetrial poolers. Their out maximum periodical daily equivance (JDEd are these in Table 1.





Thank you for your attention!

Helmut Ernstberger

Product Specialist Inorganic EMEA EMT and Spain helmut.ernstberger@perkinelmer.com

